**CENTRÍFUGA BAJA VELOCIDAD MULTIBAS+
MULTIBAS+ LOW SPEED CENTRIFUGE
CENTRIFUGEUSE BASSE VITESSE MULTIBAS+**

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Thank you for choosing this equipment. We sincerely wish that you enjoy your Nahita-Blue Multibas + centrifuge. We highly recommend looking after this equipment according to what is stated in this manual.

Nahita develops its products according to the CE marking regulations as well as emphasizing the ergonomics and security for its user. The correct using of the equipment and its good quality will permit you to enjoy this equipment for years.

The improper use of the equipment can cause accidents and electric discharges, circuit breakers, fires, damages, etc. Please read the point of Maintenance, where we expose the security notes.

TO GET THE BEST RESULTS AND A HIGHER DURATION OF THE EQUIPMENT IT IS ADVISABLE TO READ THOROUGHLY THIS MANUAL BEFORE OPERATING WITH THE EQUIPMENT.

Please bear in mind the following:

- ◆ Please read this manual carefully for the first-time using.
- ◆ Low speed centrifuge Multibas+ can only be operated by trained or authorized staff.
- ◆ The maintenance of the equipment can only be done by Nahita Technical Service Department or by authorized distributors.
- ◆ NEVER use the centrifuge with the following materials:
 - Flammable or explosive materials
 - Strengthen chemical materials
 - Toxic or radioactive substances or pathogenic micro-organisms
- ◆ Use only original accessories, the use of other type of accessories could produce damages to the equipment and the users.
- ◆ Inspections and maintenance must be conducted periodically.
- ◆ In order to avoid any injury or damage to the person, equipment or environment, please comply with all the security statements in this manual.
- ◆ Besides the measures for accident prevention, environmental protection and the security and professional career common rules, the user should obey all the rules and laws of the country and local region.
- ◆ Any action ignoring the safety procedure, the law and rules, or other relative rules will cause damage and danger to the person, device and environment.

Caution: Before using of equipment, please read the Manual carefully!

Caution: danger of low tension! Electricity!

ENVIRONMENTAL USE CONDITIONS

◆ The centrifuge may be damaged by the following factor that could compromise the safety in using the equipment:

- Chemical effects
- Environmental impacts, including radiation of natural UV
- Corrosion, abrasion and wear of shield and other security components

◆ The equipment should be placed in a location:

- Properly ventilated
- Free of vibrations and airflows that could affect the performance of the instrument
- Free of corrosive or explosive substances

◆ Indoor use;

- Altitude: ≤ 2000 m
- Temperature of use range: $+5^{\circ}\text{C} \sim +40^{\circ}\text{C}$
- Relative humidity $\leq 80\%$
- Power supply: CA 230V / 50-60 Hz / 10A

ATTENTION!! IF EQUIPMENTS ARE NOT PROPERLY CLEAN AND DISINFECTED THEY WOULD NOT BE ALLOWED TO REPAIR BY OUR TECHNICAL SERVICE.

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INDEX OF LANGUAGES

1. Security	28
1.1 Operation notes	29
2. Introduction	29
2.1 Overview	29
2.2 Outlook	29
2.3 Structure	30
2.4 Security protections	31
2.5 Placement requirements	31
3. Rotor type and technical parameters	31
4. Preliminary preparation	33
4.1 Transport and installation	33
4.2 Choose the right installation place	33
4.3 Changing location	33
4.4 Connection to power supply	34
5. Operating instructions	34
5.1 Control panel and display	34
5.2 Turn on	35
5.3 Open the lid	35
5.4 Close the lid	35
5.5 Rotor installation	35
5.6 Maintenance of rotor shaft	36
5.7 Calculation of maximum permitted rotor load	36
5.8 Loading samples	37
5.9 Safety use of rotor	37
5.10 Operation parameter setting	38
5.11 Other parameters and functions	39
6. Maintenance	40
6.1 Cleaning and purification	40
6.2 Maintenance	41
6.3 Warranty	41
7. Troubleshooting	42
7.1 Emergency lid opening	42
7.2 Information on possible failures and problems	42
8. Technical parameters	45

1. SECURITY

Multibas+ centrifuge is based on current technology and safety standards:

EN 61010-1:2010. Safety requirements for electrical equipment for measurement, control and laboratory use -- Part 1: General requirements.

EN 61010-2-020:2006. Safety requirements for electrical equipment for measurement, control, and laboratory use -- Part 2-020: Particular requirements for laboratory centrifuges.

Therefore, using must be in accordance with the requirements of the design. The incorrect or inappropriate use of the Multibas+ centrifuge will result in equipment damage or personal injury:

- Use not in accordance with design requirements
- Use by non-qualified or authorized operation and maintenance staff
- User changes content inappropriately without authorization of the design
- Not pay attention to or understand the rules of safe use

ANY PERSON INVOLVED IN USE OR MAINTENANCE OF THIS EQUIPMENT MUST READ AND UNDERSTAND THIS MANUAL AND THE SECURITY RULES SPECIFIED IN THIS MANUAL.

The following safety rules and recommendations must be taken into account to avoid accidents that could produce damages or injuries to the user of the equipment as well as damages to the equipment or sample:

- ◆ This manual is part of the Multibas+ centrifuge, so it should be available for all users of the equipment.
- ◆ The Multibas+ centrifuge has been designed to be used in biology, chemistry, genetics, immunology, etc. At a maximum separation speed the density of the sample must not exceed 1.2 g/cm³; when the sample density is higher than 1.2 g/cm³ the maximum rotor speed must be reduced accordingly.
- ◆ When the Multibas+ centrifuge is working (during the process of sample separation or with the rotor turning), the user must not stay within 30 cm around the centrifuge; also there must not be hazardous substances or items blocking vents within 30 cm around the centrifuge.
- ◆ Centrifuge design is neither anti-corrosion nor explosion-proof, therefore, ensure the environment is free of explosive or corrosive substances.
- ◆ The use of the centrifuge is strictly prohibited with the following materials:
 - Flammable or explosive materials
 - Strengthen chemical materials
 - Toxic or radioactive substances or pathogenic micro-organisms
- ◆ The separation of corrosive substances can cause damage and destruction of centrifuge internal materials or weaken the rotor mechanical strength; so in case of separation of corrosive substances they must be kept in adequate and properly protected containers.

Problem	Causes y solutions
Smell of burnt from the machine	- Immediately turn the equipment off and unplug it. - Motor or any electrical component is burnt. - Contact your distributor for repairing through Technical Service Department
Not any of the above problems	- Contact your distributor for repairing through Technical Service Department

8. TECHNICAL PARAMETERS

Environmental conditions of use

- Laboratory indoor use
- Installation place properly ventilated, free of vibrations or air flows and free of corrosive or explosive substances
- Altitude: ≤2000m
- Relative humidity: ≤80%

<i>Ambient temperature</i>	+5°C ~+40°C
<i>Power</i>	AC230V 50HZ/60HZ 10A
<i>Timer</i>	1 h 59 min 59 s
<i>Max. speed</i>	5000 rpm
<i>Max. relative centrifugal force</i>	4980 g
<i>Acceleration/Deceleration rates</i>	0-9
<i>Acceleration time</i>	< 35 s from 0 to max. speed
<i>Breaking time</i>	< 30 s from max. Speed to 0
<i>Programs</i>	9
<i>Short spin function</i>	Yes
<i>Noise (at max. Speed)</i>	≤65dB(A)
<i>Dimensions (machine)</i>	590x470x350 mm (LxAxH)
<i>Net weight</i>	32 Kg

INSTRUCTIONS ON ENVIRONMENT PROTECTION

At the end of its life cycle, please, do not dispose of this equipment by throwing it in the usual garbage; hand it over a collection point for the recycling of electrical and electronic appliances. It does not contain dangerous or toxic products for humans but a non adequate disposal would damage the environment.

The materials are recyclable as mentioned in its marking. By recycling material or by other forms of re-utilization of old appliances, you are making an important contribution to protect our environment.

Please inquire at the community administration for the authorized disposal location.

Chart 4: troubleshooting

Problem	Causes y solutions
No display	- Check whether the power outlet and the connection is right and the power supply is on - Check whether the switch is on - Check the fuse - If problem still persist, contact your distributor for repairing through Technical Service Department
Stop to work suddenly	- The working speed exceeds the maximum allowed speed of rotor - If working speed exceeds in 450 rpm the maximum allowed speed of rotor, the alarm of over speed will sound. Wait until the centrifuge stops and reset the speed - Working speed exceeds selected speed - If motor is overheated, power supply will be cut inside the equipment and then it will stop - If display is not on, check the power supply - Maybe input voltage is too low; please check whether the voltage of the power supply is right
Lid cannot be open	- Lid cannot be open until rotor completely stops - Open the lid manually - If problem persists, contact your distributor for repairing through Technical Service Department
Strong vibration of equipment	- This situation is normal when rotor reaches the critical speed - Check whether the rotor is locked tightly or not - Check the symmetry of the load of the rotor; check the working situation of the machine - Check the motor shaft: if it cannot be rotated manually and smoothly then maybe there is a problem with shaft of motor; contact your distributor for repairing through Technical Service Department
Front panel is not on after the equipment is switched on	- Power supply is not connected, please check the power - Fuse is blown, replace the fuse - If problem persists, contact your distributor for repairing through Technical Service Department
Unusual display of operating panel	- Maybe it's because of the interference of power grid, please shut down and start it after 1 minute, and then the display will be right. - If problem persists, contact your distributor for repairing through Technical Service Department.
Motor doesn't work	- Electric control circuit is broken; contact your distributor for repairing through Technical Service Department

1.1 OPERATION NOTES

- ◆ Before centrifugation, make sure the rotor used is adequate and it is correctly installed.
- ◆ During centrifugation (rotor running) or during stop process (but with rotor still rotating), do not open manually the lid or move the equipment.
- ◆ Use only the original accessories supplied with the equipment. For the use of other common accessories as glass or plastic containers, it is necessary to make sure they meet with the requirements of use as resistance to speed or RCF used, etc.
- ◆ Never open the lid during operation.
- ◆ The repair or replacement of mechanical parts and electronic devices of the centrifuge, must be done only by Nahita Technical Service Department.
- ◆ For a correct use of the equipment it is necessary to select the appropriate rotor load as well as avoid overloading.
- ◆ Ensure regular checks of the rotor; if rotor has obvious corrosion or damage signs it must be replaced immediately.
- ◆ Maintenance should be done after a period of use, according to point "Maintenance".


2. INTRODUCTION

2.1 OVERVIEW

Multibas+ centrifuge has been designed to be used in biology, chemistry, genetics, immunology and other fields in which the use of conventional centrifuges for separation and sedimentation of samples is required. For configuration of the equipment please refer to Point 3 "Rotor type and technical parameters".


2.2 OUTLOOK

figure 1: front side

1. Lid
2. Shell
3. Legs
4. On/Off switch
5. Emergency opening system
6. Front case
7. Air intake
8. Control panel and display
9. Peephole
10. Power socket

figure 2: back side

2.3 STRUCTURE

This equipment consists of the following parts: cover, shell and base, inner chamber, mechanical components (motor and rotor), electrical components (power, display and control systems) and other components as alarm system.

1. The cover system includes the lid, hinges and springs, locking system, alarm and emergency opening. The hinges are at the back side and the locking system at the front side; only when the lid is closed the centrifuge can start running, otherwise the alarm system (sound) will activate and operation will not be able to start.

◆ To open the lid, press the lid open button  on the control panel. When the lid is lifted manually up to a certain height, the hinges and springs are able to support its weight.

◆ In case of power failure or problems with the open button  and if it is necessary to take out the samples on time, it will be necessary to use the emergency opening system; for this, pull the cord slowly until the lid is open (the position of emergency lock cord is show non *figure 1*).

Note: never open the lid manually with the equipment running or turned on!

2. The inner chamber consists of a stainless steel chamber and sealing rubber joints. This inner chamber provides a stable working environment.

3. The equipment uses a variable frequency electric motor to move the rotor that contains the samples. The motor shaft connects with the rotor by means of a conical system that assures a smooth and precise movement.

4. The rotor system includes adapters for tubes of different capacities (for more detail, please refer to Chart 1. "Rotors, type of adapters and technical parameters") and other related accessories. The function of rotor is to serve as support of samples for centrifugation at a certain rotation speed that results in a certain relative centrifugal force that produces the separation of sample. As a result of the rotation speed, the relative centrifugal force is thousands of times the Earth gravitational acceleration (g), so it is essential a safety use and maintenance of the rotor.

5. The base of the equipment is composed of a steel protecting plate and supporting legs.

6. The power system includes the power socket and switch, which is responsible of the equipment's electrical supply from the power net.

7. The control system includes setting of speed, relative centrifugal force, time, acceleration/deceleration rates and the whole system of display, alarms, etc. To assure a proper and safety use never disassemble the equipment.

8. The display system consists of a LCD display and PVC keyboard (control panel). This is the interface for user-machine dialogue. It can simultaneously display and set the different parameters as well as display a variety of errors.

9. The alarm systems includes lid cover fault and over speed. When the equipment detects a fault or problem, the alarm system will start with buzzer sound, the red light on display will light and centrifuge cannot start running or it will stop automatically; once problem is solved, the centrifuge can be re-started. Note: to mute the alarm sound, press the start/stop button  on the control panel.

Chart 3: Error symbols and their causes

Symbol	Cause	Solution
Erro1	Imbalance: centrifuge stops due to a detection of an excess of vibrations	- Re-balance the samples; the maximum error allowed is $\leq 1.5g$ - The device is placed in a leaning position, not horizontal; change position of equipment so as it is perfectly horizontal. - The motor shaft is bent; equipment needs to be repaired by Technical Service Department
Erro2	Over speed: centrifuge stops due to a detection of over speed of motor	- Problem of microcomputer control system; equipment needs to be repair by Technical Service Department. - Problem of speed sensor; equipment needs to be repaired by Technical Service Department
Erro3	Lid protection system is acting	- Re-close the lid - Lid protection system is broken; equipment needs to be repaired by Technical Service Department
Erro4	Input voltaje too low	- Check the power source
Erro5	Braking over voltage	- The Power source voltage is too low or the braking resistor failed or brake too fast; equipment needs to be repaired by Technical Service Department
Erro6	Over current	- The Power source voltage is too low or the acceleration is too fast; equipment needs to be repaired by Technical Service Department
Erro7	No speed measurement	- Speed sensor or motor failed; equipment needs to be repaired by Technical Service Department
Erro8	Fault of communication between internal components of the equipment.	- Equipment needs to be repaired by Technical Service Department
Erro9	Over voltage	- The Power source voltage is too low or the acceleration is too fast; equipment needs to be repaired by Technical Service Department

- ◆ Neither fuses nor accessories (including their loss), are covered by the product's warranty. The warranty neither covers piece's deterioration due to the course of time.
- ◆ Please, make sure you keep the invoice either for having the right to claim or asking for warranty coverage. In case you have to send the equipment to Technical Service Department you should enclose the original invoice or a copy as guarantee.

7. TROUBLESHOOTING

7.1 EMERGENCY LID OPENING

- ◆ In case of accidental power failure the lid will not be able to open automatically; if it is necessary to remove samples from the inner chamber, the lid should be open manually.

Note: this opening system is only for emergencies; do not use it as usual way.

- ◆ If a power failure happens, the motor will take longer time to stop completely since the braking function will be also useless. Please, wait patiently until the motor completely stops.

For emergency opening, proceed as follows:

- ◆ Confirm the motor has completely stopped.
- ◆ Unplug the equipment.
- ◆ Pull firmly and horizontally from the emergency cord; the lid will open and samples could be extracted from the centrifuge.

7.2 INFORMATION ON POSSIBLE FAILURES AND PROBLEMS

The below chart shows possible failures and their causes and solutions. In spite of the indicated recommendations the problem still persists, please contact your distributor for the reparation through the Nahita Technical Service Department.

PLEASE, IN CASE OF ANY PROBLEM TURN THE CENTRIFUGE OFF AND DISCONNECT IT FROM THE NET.

2.4 SECURITY PROTECTIONS

Multibas+ centrifuge has a series of security protections:

- ◆ The base is made of steel and the inner chamber is made of stainless steel.
- ◆ The lid uses explosion-proof structure, has a peephole in the central part and has locking system. Only when the centrifuge is turned on and the rotor is stopped, the open button  can be pressed to open the lid; only when the lid is locked, the centrifuge can start running.
- ◆ Over speed: when the rotor speed exceeds in 250 rpm the set speed, the alarm will sound; when the rotor speed exceeds in 300 rpm the maximum permitted rotor speed, the motor will automatically stop. Once rotor is completely stopped, open the lid, solve the problem and re-run.
- ◆ Emergency lid opening system: if during operation there is a sudden power failure or machine fault and the lid cannot be open through the button on control panel then it will be necessary to use the emergency opening system (please refer to "Troubleshooting").

2.5 PLACEMENT REQUIREMENTS

- ◆ The centrifuge must be placed on a rigid, stable surface and must be kept away from any impact or vibration; avoid heat and direct sunlight exposure.
- ◆ There must be a free space of 10-15 cm all around the centrifuge to assure correct cooling ventilation.
- ◆ Level must be adjusted after installation and make the four supporting legs at the bottom of the equipment to support evenly on the horizontal surface.
- ◆ Power supply to the equipment must be of 230V AC / 50-60HZ and must be properly earthed.

Note: the equipment must be perfectly grounded. Do not cut off power supply during operation, the control circuit could be damaged.

3. ROTOR TYPE AND TECHNICAL PARAMETERS

Multibas+ centrifuge is supplied with swing-out rotor and can be used with other swing-out rotors and different adapters for different type of tubes.


figure 3: swing-out rotor and different type of adapters


figure 4: rotor for microplate

Chart 1: Rotors, type of adapters and technical parameters

Rotor	Rotor reference	Adapter reference	Code*	Capacity (mLxn° tubes)	Max. speed (rpm)	Max. RCF (xg)	Type of tube
Swing-out (included) • (adapters not included)	—	92751002	#5001	50×4	5000	4980	PP conical bottom (Falcon type)
		92751001	#5002	100×4		4600	PP
		92751003	#5003	50×8	4000	3040	PP conical bottom (Falcon type)
		92751005		15×32			
		92751004		15×24			
		92751007		7×48			
		92751008		7×72			
		—	#5004	2×2×96 wells	3500	2120	PP
		—	#5005	250×4	5000	4700	PP
Swing-out for microplates (optional)	92751009	—					
Swing-out (optional)	92751000	—					

*Code to be selected on the control panel of the centrifuge

6.2 MAINTENANCE

- ◆ During installation always handle the rotor carefully avoiding hits or knocks; small damages in the rotor could cause breakage during centrifugation.
- ◆ Check the rotor and adapters regularly to detect possible defects, damages, breaks or knocks; if any damage is detected replace the component immediately.
- ◆ To remove the rotor, hold it firmly and lift it vertically. Avoid sudden movements.
- ◆ Under normal circumstances the rotor should be cleaned once a week; in case of using salts or corrosive substances containing samples, cleaning should be made after each use as well as in case of spillages or splashes of the centrifuged solutions.
- ◆ For cleaning the rotor use neutral detergents or a sponge or soft cloth and then wipe with distilled water. Do not spray water on the rotor since liquid could stay or reach other components and cause corrosion.
- ◆ Use a cloth or tweezers to remove the solid residues from the inner chamber of the centrifuge.
- ◆ In case of breakdown do not try to repair the equipment by yourself; apart from losing warranty you can cause damages in the general performance and operation of the equipment as well as damages to people (burns, wounds,...) and the power supply or nearby electrical devices.
- ◆ Only use original accessories, the use of other type of accessories could cause damages to the equipment and the users.
- ◆ Make sure the centrifuge is turned off and disconnected from the net in case of long periods of no use.

Maintenance of the motor shaft

- ◆ Turn the centrifuge on through the general switch and wait until self-test is finished.
- ◆ Press button  to open the lid.
- ◆ Use the supplied wrench to disassemble the rotor and remove it from the motor shaft; clean the rotor nut. Note: turn clockwise to tighten the nut and anti-clockwise to loosen it.
- ◆ Clean the conical part of the motor shaft, and then grease it well.

Transport, storage

The centrifuge is a precision equipment so during transport or storage prevent it from moisture, knocks or horizontal positions.

6.3 WARRANTY

- ◆ This equipment is protected under the Warranties and consumer goods regulation (10/2003).
- ◆ Overhaul is not covered by the equipment warranty.
- ◆ Operations made by non-qualified staff will automatically produce a loss of the equipment warranty.

5. Calculation of relative centrifugal force (RCF)

Relative Centrifugal Force (RCF) is the necessary force to produce separation of particles during centrifugation. The unit of this force is expressed in number of times the gravity (xg) and its value is directly related to centrifugation speed and centrifugal radius.

The RCF is calculated by the following formula:

$$F.C.R. = 1.118 \cdot 10^{-5} \cdot r \cdot n^2$$

Being $1.118 \cdot 10^{-5}$ a constant.

r = centrifugal radius, horizontal distance in centimeters from the rotating axle up to the bottom of the tube (cm).

n = rotating speed expressed in revolutions per minute (rpm).

6. MAINTENANCE

6.1 CLEANING AND PURIFICATION

- ◆ In case of spills of hazardous substances inside the equipment, user has the responsibility to make an appropriate cleaning and purification.
- ◆ User should perform a correct cleaning and purification of equipment according to what is stated in this manual to avoid damages to the centrifuge. The use of inappropriate products and/or cleaning procedures can cause damages to the equipment and in its internal components.
- ◆ Turn the equipment off and disconnect it from the net before cleaning.
- ◆ A cleaning of the case, inner chamber, rotor and adapters must be done regularly to avoid accumulation of residues that cause damage or corrosion of the different components.
- ◆ Do not use organic solvents that could remove the lubrication of the gears and bearings of the motor.

Steam sterilization

In *Chart 2*, the life time period of different centrifugation tubes is shown in relation to the number of sterilization cycles and use. If any damage is detected on rotor, adapters or centrifugation tubes, replace immediately.

	Max. temperature (°C)	Min. period of time (min)	Max. period of time (min)	Max. number of cycles
Tube glass	134 - 138	3	5	—
Tube PC	115 - 118	30	40	20
Tube PP	115 - 118	30	40	30
Tube PA	115 - 118	30	40	20

Chart 2: sterilization parameters

4. PRELIMINARY PREPARATION

4.1 TRANSPORT AND INSTALLATION

- ◆ Multibas+ centrifuge is supplied in a box with all the necessary protections to assure a safety transport of the equipment. Please, open the box and remove all protections.

The net weight of the equipment is approximately 30 Kg; to avoid damages in the motor shaft, remove the rotor before moving or transporting the centrifuge. Do not shake the centrifuge.

- ◆ For long-distance transport or movement, please use the original box, close it properly and keep it always in vertical position; handle the equipment with care.

4.2 CHOOSE THE RIGHT INSTALLATION PLACE

- ◆ Multibas+ centrifuge has been designed for indoor use and the installation place should be in accordance with the below requirements.
- ◆ When the centrifuge is running, a safety distance of 30 cm around the equipment must be kept free of hazardous substances and free of people.
- ◆ The supporting surface must be solid, horizontal and free of movements or vibrations; in case of using movable surfaces, the centrifuge must be perfectly held to assure a safety operation.
- ◆ There must be a minimum distance of 10-15 cm between back and lateral sides of the centrifuge and the wall or adjacent equipments in order to ensure smooth flow of air and cooling cycle of the equipment.
- ◆ Centrifuge should be placed away from windows to avoid heat and direct sunlight exposure.
- ◆ Once the equipment is placed, make sure the four legs stand evenly on the horizontal surface and the equipment is well level.
- ◆ The room in which the centrifuge is placed must be clean and with a constant room temperature between 5-40°C and a relative humidity ≤ 80%.

4.3 CHANGING LOCATION

Once the centrifuge has been installed, it should not be moved. In case of being necessary to change location, the new location should be solid, horizontal, free of movements and vibration and the centrifuge should be perfectly level.

To avoid damages on the rotor shaft, please, remove the rotor before moving the equipment.


4.4 CONNECTION TO POWER SUPPLY

♦ The power supply connected to the centrifuge must be properly grounded. Check the cable used meet with the safety requirements in your country or region. The inlet voltage and frequency should match with the requirements indicated in the plate of specifications of the centrifuge (230V / 50-60Hz / 10A). Please, use the power cable supplied with the equipment to connect the centrifuge to the power net.

♦ The equipment will be turned on with the switch on position "I" and turned off with the switch on position "O".

5. OPERATING INSTRUCTIONS

5.1 CONTROL PANEL AND DISPLAY

- | | |
|--------------------------------------|---------------------------|
| 1. Lid open button | 11. Rotor/adapters code |
| 2. Stop button | 12. Set acceleration rate |
| 3. Start/Short spin button | 13. Set deceleration rate |
| 4. Parameter setting knob | 14. Program number |
| 5. Temperature inside the centrifuge | 15. Lid status |
| 6. Software version | 16. Error code |
| 7. Remaining operation time | 17. Set operation RCF |
| 8. Set operation time | 18. Real operation RCF |
| 9. Real operation speed | |
| 10. Set operation speed | |


figure 7

6. In case of mechanical failure or not correct parameter during parameter setting process, an alarm will sound. Press button to stop the alarm and review the selected parameters.

7. Press button , the centrifuge will start running, time will gradually decrease and when zero is reached the equipment will automatically stop (if the centrifuge needs to be stopped during on half-way press the button). Speed will also decrease (related to the selected deceleration rate) and when zero is reached the centrifuge lid will automatically open; an alarm of end of centrifugation will sound and it can be stopped by pressing . The centrifugation is completed.

5.11 OTHER PARAMETERS AND FUNCTIONS

1. **Centrifugal force setting:** to set RCF instead of speed press vertically the parameter setting knob until digits corresponding to set RCF blink. Then turn the parameter setting knob until the desired value is displayed.

2. **Programs:** the centrifuge has the possibility of saving up to 9 programs. During parameter setting, once all operation parameters have been selected, select the program number and automatically parameters will be saved in that program number.

3. For short spin of samples keep pressed button , and the centrifuge will start running until reached the selected speed and centrifugation time will be displayed. Once button is released the centrifuge will stop.

4. Once speed becomes stable, the different parameters of RCF, speed, time or acceleration and deceleration rates can be modified. After modification it is not necessary to confirm manually since digits will flash 3 times to confirm they have been saved.

Note: before configuration of operation parameters it is necessary to install correctly the rotor on the motor shaft. In cases of errors in parameter setting process, it is possible to re-set the parameters.

DO NOT START USING THE EQUIPMENT UNTIL THE NUT IS PERFECTLY TIGHTEN TO THE MOTOR SHAFT.

- ◆ After a continuous use of the equipment, check if the rotor nut is loosen and tighten it if necessary before next use.
- ◆ Sample tubes must be placed in rotor at the same time and in symmetrical positions (permitted difference is $\leq 1.5\text{g}$); asymmetrical sample load is not permitted.

5.10 OPERATION PARAMETER SETTING

Taking as example the use of swing-out rotor with adapters 8x50 mL (#5003 according to *Chart 1*), proceed as follows for parameter configuration. For example, value of parameters will be the following:

Code	Speed (rpm)	Time (min/s)	Acceleration	Deceleration
#5003	4000	30 min 00 s	5	3

1. Firstly turn the equipment on by pressing the On/Off switch placed at one side of the centrifuge; the LCD display will light on.
2. **Speed setting:** Press vertically the parameter setting knob. Digit corresponding to set Speed will blink. Turn the parameter setting knob until number 4000 is displayed.
3. **Time setting:** Press vertically the parameter setting knob. Digit corresponding to set Time will blink. Turn the parameter setting knob until 30 min 00 s is displayed.
4. **Rotor/adaptor number setting:** Press vertically the parameter setting knob. Digit corresponding to Rotor/adapters code will blink. Turn the parameter setting knob until #5003 is displayed. Once operation parameter has been set press the setting knob vertically or wait until value on display flashes 3 times indicating that the value has been saved.
5. **Acceleration and deceleration rates setting** (value from 0 to 9), increasing values reduce the time for reaching setting speed or zero speed respectively). Press vertically the parameter setting knob. Digit corresponding to set Acceleration rate will blink. Turn the parameter setting knob until number 5 is displayed. Press again vertically the parameter setting knob and digit corresponding to set Deceleration rate will blink. Turn the parameter setting knob until number 3 is displayed.

Once operation parameter setting procedure is finished, display will be as *figure 7*:

5.2 TURN ON

- ◆ Connect one end of the cable supplied in the socket of the centrifuge placed at the back side of the equipment and the other end to the power net (AC230V 50HZ/60HZ) properly grounded.
- ◆ After connecting the equipment, press the switch at the back side of the equipment. With a short buzz, the LCD display on control panel will turn on and after a complete auto-test the equipment will be ready to go to next point.

5.3 OPEN THE LID

- ◆ Press the lid open button  placed on the control panel and the lid will be automatically open; lift the lid manually up to a certain height where it can be held by the spring and hinges and making accessible the inner chamber. If it is the first time opening the centrifuge remove the rotor and its packaging.

Note: When pressing the lid open button, make sure not to have the head or other objects above the lid to avoid hits.

- ◆ In case of failure and the button does not work, it will be necessary to open the lid manually according to point "Troubleshooting".

5.4 CLOSE THE LID

- ◆ Push the lid downwards until the hook at the front of the lid passes the lock pin with a "click" sound; then the bottom of hook will touch the trip switch and the lid will be locked.
- ◆ Push the lid downwards properly; do not overexert or the hook will be damaged.

5.5 ROTOR INSTALLATION

- ◆ It is only allowed to use original rotors and adapters designed and manufactured for Multibas+ centrifuge, otherwise the equipment could be damaged; the rotors and adapters allowed are specified in this manual (refer to Point 3 *Rotors, type of adapters and technical parameters*).

THE USE OF INAPPROPRIATE ROTORS, ADAPTERS OR TUBES MAY DAMAGE THE CENTRIFUGE.

◆ Proceed as follows for rotor installation (see figure 6):

◆ Turn the equipment on and wait until auto-check is finished.

◆ Press the open button  to open the lid; check the inner chamber is clean and free of any object.

◆ Clean the surface of the motor shaft.

◆ Pass the central hole of the rotor vertically through the conical surface of the motor shaft and press until it is perfectly adjusted. The conical surface allows a wear-free adjustment and rotation of rotor on the motor shaft.

◆ By using the hexagonal wrench tighten firmly the rotor nut on the motor shaft and make sure it does not move; otherwise adjust it correctly.


figure 6

Note: once rotor is installed, each time before and after use it is necessary to check if rotor position has changed or not and tighten the nut if necessary to assure a solid installation of rotor.

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5.7 CALCULATION OF MAXIMUM PERMITTED ROTOR LOAD

There is a big centrifugal force when Multibas+ centrifuge is working at maximum speed. The rotor is designed with the enough mechanical strength to support the maximum permitted speed; this is called "security factor". However, this security factor requires that the rotor load does not exceed the maximum permitted load.

If the sum of samples, tubes and adapters exceed the maximum permitted load of rotor then it is necessary to reduce the weight of the sample or re-calculate the rotor permitted speed (NPERM) to assure that the rotor load does not exceed the maximum permitted load:

$$NPERM = N_{max} \times (\text{Max. Permitted load/real load}) \times 0.5$$

For this, it is necessary to take into account that the density of sample must not be higher than 1.2 g/cm³; if density of sample is higher than this value then the maximum permitted speed of rotor should be reduced and re-calculated according to above formula.

DON'T USE THE ROTOR OVERLOAD, OR IT WILL CAUSE EXPLOSION OF THE ROTOR AND THE DEBRIS WILL DAMAGE THE CENTRIFUGE

5.8 LOADING SAMPLES

If samples are well balance when loaded in the centrifuge, separation of the different phases during centrifugation will be optimal due to the absence of vibrations that could interfere. Therefore, the load of samples must be as much balanced as possible by placing tubes of same weight, shape and size in geometrically symmetrical positions; if necessary use tubes with water. This will avoid possible vibrations that will assure a better performance in separation of samples.

The selection of adapters and tubes for centrifugation must be done taking into account the maximum nominal value of speed or RCF (see Chart 1); do not exceed these values.

Take into account the using life time of the recipient for simple centrifugation (plastic, glass). Check if they present any damage and if affirmative replace by a new one.

5.9 SAFETY USE OF ROTOR

◆ Sample loading and tube placing must be perfectly balanced and symmetrical.

◆ Once adapters are placed in the corresponding position of the rotor, check by hand that they swing freely.

◆ Swing-out rotor must not run at the critical speed of 800 rpm for a long time or equipment will suffer big vibrations that will affect its service life.

◆ In case of rotor replacement, use the supplied wrench to loosen the rotor nut by anti-clockwise rotation and then replace rotor.