MANUAL DE INSTRUCCIONES

AGITADOR ORBITAL O VAIVÉN

RSLAB-7 (2,5kg) agitador orbital, ref. LJN003

RSLAB-7 (2,5kg) agitador vaivén, ref. LJQ002

RSLAB-7PRO (7,5kg) agitador orbital, ref. LJN004

RSLAB-7PRO (7,5kg) agitador vaivén, ref. LJQ003

ÍNDICE

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




500ml 9

INSTRUCTIONS MANUAL

LINEAL AND ORBITAL DIGITAL SHAKERS

RSLAB-7 (2,5kg) orbital shaker, code LJN003

RSLAB-7 (2,5kg) lineal shaker, code LJQ002

RSLAB-7PRO (7,5kg) orbital shaker, code LJN004

RSLAB-7PRO (7,5kg) lineal shaker, code LJQ003

INDEX

1. Safety Instructions	p.3
2. Use	p.4
3. Inspetion	p.4
4. Trial Run	p.4
5. Control and Set	p.4
a. Control	p.4
b. Set Time	P.5
c. Setting Speed	p.5
6. Operating Modes	p.5
a. Mode A	p.5
b. Mode B	p.6
c. Switching the Mode	p.6
7. Supported load	p.7
8. Trouble shooting	p.7
9. Maintenance and Cleaning	p.7
10. Normas	p.7
11. Technical Data	p.8
12. Accessories	p.8

This is the User Manual of the Lineal and Orbital Digital Shakers RSLAB-7 and 7PRO. In case of need, please contact your distributor and he will help you and advise for the correct use of the equipment. It will be required that the user provide the serial number of the equipment, describe the problems, the procedures made in order to solve the problems and the person contact data.

RSLAB stirrers have a 12 months warranty from the date of the invoice for manufacture faults in case of a correct use as described in this manual. This warranty is not applied to the stirrer or any spare part as consequence of a wrong installation, wrong connections, incorrect use, accident o conditions not according to the proper use.

Contact your distributor for any complaint under warranty. The shipment only is possible by previous agreement with the distributor. In this case, it is required attached the mail explaining the problems detected.

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Warning !

Read the operating instructions completely before use and follow the safety instructions.
Ensure that only trained staff work with the instrument.

Protective ground contact!

Make sure the socket must be earthed (protective ground contact) before use.

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Use the instrument following the safety instructions and profession safety rules to avoid any accident.

Users are not allowed to touch moving parts of he instrument to avoid any accidents. Pay attention to protect hands or fingers from rolling when moving the instrument.

Set up the instrument in a spacious area on an stable, clean, non-slip, dry and fireproof surface. Do not operate the instrument in explosive atmospheres, with hazardous substances or under water.

Observe the vessel on the shaking table for medium splashing out of liquid when setting the motor speed. Reduce the motor speed if the instrument is not running smoothly.

No utilice componentes dañados. Safe operation is only guaranteed with the accessories described in the "Accessories" chapter and by following the operating instructions in full. Always disconnect the plug before fitting accessories. Check the instrument and accessories beforehand for damage each time you use them. Do not use damaged components

Do not cover the instrument. Protect the instrument and accessories from bumps and impacts. Keep away from high magnetic field.

The operating voltage required as stated on the instrument and voltage in the supply network must match.

The instrument may only be opened by trained technicians. And always disconnect the plug beforehand.

2- Use

The instrument is designed for mixing or shaking liquids in schools, laboratories or factories. Shakers RSLAB-7 y RSLAB-7PRO is suitable for usage in various attachments for mixing liquids in bottles, flasks,

test tubes for a maximum supported weight of 7.5 kg (330) and 2.5Kg (180). This instrument is not suitable for use in residential areas or other areas that may cause danger to the user or instrument as mentioned in chapter 1.

3- Inspection

Unpack the instrument carefully and check for any damages which may have arisen during transport. If it happens, please contact distributor for technical support.

If there is any apparent damage to the system, please do not plug it into the power line.

RSLAB-7 and RSLAB-7PRO are supplied with the following items:

-Stirrer	1ud
-Power cable	1ud
-User manual	1 ud

4- Trial run

Make sure the required operating voltage and power supply voltage match.

Ensure the socket is securely earthed.

Plug in the power cable and power ON the instrument.

LCD displays the safe rotary speed limit and operating mode.

Turn the speed knob to set the rated rotary speed


Press the speed knob. And the shaking function is switched ON.

Press the motor regulation knob again. and the shaking function is switched OFF.

If these operations above are normal, the instrument is ready to operate following the operation" chapter.


If these operations are not normal, the instrument may be in security protection state or be damaged. Please contact with distributor Technical support.

You can improve the stability of the instrument with the adjustable feet.

5- Control and set

a. Control

Items	Decription
Speed Knob	Set the rated rotary speed. "rotation" is switched ON or OFF via pressing the knob.
Time Knob	Set working time. "rotation" is switched ON or OFF via pressing the knob
Shaking table	Plato de agitación para los recipientes
LCD Speed / LCD Time	Set working time. "rotation" is switched ON or OFF via pressing the knob

Put the instrument on stable and safe place and plug in the mains power.

Turn ON the mains switch on the left panel. The instrument starts self-test. The instrument shows rated speed and time after initialization..

Turn the speed knob on the left side to set the rated speed. Press the speed knob and characters on the left LCD do not flash any more, and the shaking function is switched ON. Press the speed knob again and characters on the left LCD begin to flash, and the shaking function is switched OFF.

Turn the time knob on the left side to set the rated time. Press the time knob and characters on the LCD do not flash any more, and the timing function is switched ON. Press the time knob again and characters on the right LCD begin to flash, and the timing function is switched OFF..

b. Set time

The equipment can work in continuous mode or time/continuous (max. 19h59m). The setting of the work time is made through the control in the right side. In case the equipment turns off and on again, the setting disappears. During the process of running at the time and speed set up, the user can activate these two functions through the regulation controls, as indicate below:

- Pressing the time knob, the speed and the time will stop. Pressing the time knob again, restart the shake function, the time will use the pre-set time.
- Pressing the speed knob, the speed and the time will stop; Pressing the time knob again, restart the shake function, the time not(There are flashing on the right of the LCD), pressing the time knob again restart the time function. The time will use the pre-set time..

c. Setting speed

User can set and increase the shaking speed by turning the regulation control.

The continuous mode is activated pressing the control on the left without time setting. The shaking stops pressing again this control.

Increasing the shaking speed by turning the control progressively in order to avoid shocks in the equipment.

Shaking speed can be modified at any momento.

6- Operating modes

a. Mode A

Operating mode A is the default mode set in factory.

After turning ON the instrument with the mains switch, the shaking and timing functions are turned OFF. The LCD displays the set time and speed values. They are adopted or varied when the corresponding functions are turned ON.

After turning ON the instrument, the following process will be shown on the display:

The left LCD displays "SAF" and the right displays the upper speed limit (rpm) which can be set by pressing the speed knob and turning at the same time.

The LCD displays run modes "StA", "A" or "B" for about 2 seconds.

The set or stored speed value is shown in the left display, and the set or stored time value is shown in the right display. Turn the speed knob and time knob can set the speed and time.

After pressing the speed knob or time knob, the instrument begins to run at the set speed. The actual speed and time remaining are shown in the LCD. If the remaining has reached zero, the shaker stops its motion. Pressing the speed knob alone, the instrument will work at continuous mode (The time function is OFF)

b. Mode B

When the instrument is turned ON and values for speed and time are erased to zero, users can set.. The upper speed limit that has been set is taken over from operating mode A and cannot be modified.. After the power has been turned OFF, the instrument will no longer automatically start up in operating Mode B by itself.

After turning ON the instrument and selected run mode, the following process will be shown on the display.


The left LCD displays "SAF" and the right displays the upper speed limit (rpm) which cannot be changed.

The LCD displays run modes "StA" for about 2 seconds.

The set speed value is shown in the left display, and the set time value is shown in the right display. Then set the motor rotary speed and time.

After pressing the speed knob or time knob, the instrument begins to run at the set speed. The actual speed and time remaining are shown in the LCD. If the remaining has reached zero, the shaker stops its motion.

c. Switching the mode

Hold down both knobs and turn ON the instrument with the mains switch. After about 5 second, you can let go of the knobs..

7- Supported load

In order to ensure safe, the shaker must only be operated within the supported: 2,5kg or 7,5kg. Make sure the holding surface of the shaker is always clean and level. Ensure the individual shaking vessels must be placed evenly and fastened securely in the middle of the shaking table no matter one or multiple shaking vessels are placed.

8- Trouble shooting

- The instrument doesn't work:
- Check whether the power cable fitted well
 - Check whether the fuse is broken
- The speed cannot reach set value:
- Check whether it is overloaded.
- The motor does not start via pressing the speed knob and time knob:
- Check whether the time is set to zero.

9- Maintenance and cleaning

Use and maintain the instrument in right way to keep it in good operating state in order to lengthen its life time and avoid the damage to the instrument. Only use the cleaning agents which have been approved by manufacturer to clean the instrument.

Only use cleanser that we advised as below:

- Dyes: Isopropyl alcohol
- Construction materials: Water containing tenside / isopropyl alcohol
- Cosmetics: Water containing tenside /isopropyl alcohol
- Foodstuffs: Water containing tenside
- Fuels: Water containing tenside

For materials which are not listed, please request information from your supplier. Before using another than the recommended method for cleaning or decontamination, the user must ascertain that this method does not destroy the instrument. Wear the proper protective gloves during cleaning of the instrument.

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Electrical instrument may not be placed in the cleansing agent for the purpose of cleaning. The instrument must be cleaned and put it into the initial packaging carton before sending to service for repair, avoiding the contamination of hazardous.

The equipment should be clean before sending to the distributor for repair.

Switch OFF the instrument and put it in the dry, clean, stable place at room temperature in long-term disuse.

10- Regulations

Construction in accordance with the following safety standards:

EN61010-1

UL3101-1

CAN/CAS C22.2(1010-1)

EN 61010-2-10







Construction in accordance with the following EMC standards:

EN61326-1

11- Technical data

Voltage	100-240 VAC	Speed display	LCD display
Frequency	50/60 Hz	Timer	Yes
Power	30W	Timer display	LCD display
Shaking movement	Orbital/ linear	Time setting range	1-1199min
Diametro orbital	10mm (RSLAB-7PRO) 4,5mm (RSLAB-7)	Run type	Time / Continuous operation
Max. shaking weight (with accessories)	7,5kg (RSLAB-7PRO) 2,5kg (RSLAB-7)	Dimensions	420x370x100mm (RSLAB-7PRO) 340x300x100mm (RSLAB-7)
Motor type		Weight	13,5kg (RSLAB-7PRO) 8,1kg (RSLAB-7)
Motor rating input	28W	Ambient temperature	5-40°C
Motor rating output	15W	Relative humidity	80%
Permissible ON time	100%	Protection class acc. To DIN EN60529	IP21
Speed range	100-350rpm (50685180) 100-500rpm (50685080) 100-350rpm (50685130) 100-800rpm (50685030)	RS232 interface	Yes

12- Accessories

	Code	Description
	LJS001	Plate with 4 bars (RSLAB-7PRO)
	LJS011	Fixing clip attachment (RSLAB-7PRO)
	LJS012	Anti-slip platform with 8 rubbers (RSLAB-7PRO)
	LJS008	Plate with 3 bars (RSLAB-7)
	LJS009	Fixing clip attachment (RSLAB-7)
	LJS014	Anti-slip platform with 3+3 rubbers (RSLAB-7)
	LJS001	Fixing clip for Erlenmeyer 25 mL

	LJS002	Fixing clip for Erlenmeyer 50 mL
	LJS003	Fixing clip for Erlenmeyer 100 mL
	LJS004	Fixing clip for Erlenmeyer 200/250 mL
	LJS005	Fixing clip for Erlenmeyer 500 mL
	LJS006	Universal bar for LJS008
	LJS007	Universal bar for LJS010


Number of fixing clips possible per shaker:

	25ml	50ml	100ml	250ml	500ml
RSLAB-7	36	18	18	5/8*	5/8*
RSLAB-7PRO	56	28	28	13	9


*Using only 3 fixing points and no 4

Position of the fixing clips in the equipment RSLAB-7:


25ml 36


50ml 18


100ml 18


200ml 8 (Three fixed allowed)


200ml 5 (Three fixed not allowed)

500ml 8 (Three fixed allowed)


500ml 5 (Three fixed not allowed)





Position of the fixing clips in the equipment RSLAB-7PRO:

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