**CONCENTRADOR DE OXÍGENO AXAVET OLV-5
AXAVET OLV-5 OXYGEN CONCENTRATOR
CONCENTRATEUR D'OXYGÈNE AXAVET OLV-5**

Ref. ZMC017

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10. TRANSPORTATION AND STORAGE

10.1 Transportation and storage precautions

- Before transport or storage, discard the water in the humidifier bottle.
 - During transport and handling, the equipment must be kept in an upright position; inverting or laying in a horizontal position is not permitted.
 - When the storage temperature is below 10 °C, place the concentrator in a normal working environment for 8 hours before use.
- If the equipment has been stopped for a long time, it should be turned on for inspection before using it again to confirm that all functions are normal.

10.2 Requirements for the storage and transportation environment

- Ambient temperature: -20 °C ~ 50 °C
- Relative humidity: ≤ 95%
- Atmospheric pressure: 500hPa~1060hPa

10.3 Transportation

- The fully packed oxygen concentrator shall be protected against violent collisions and direct contact with rain or snow during transportation.

10.4 Storage

- Should be stored in well-ventilated indoor locations, out of intense sunlight and corrosive gases.

11. ACCESSORIES

- Do not use parts and accessories not recognised by the manufacturer to avoid adverse effects on the safety and performance of the product.
- The oxygen concentrator must use a dedicated humidifier bottle; to replace it, contact your dealer to purchase one.
- The nasal cannula supplied is a disposable product, for one patient only. The user may purchase qualified cannulae with the medical device registration certificate, referencing the model and specification of the cannula supplied.

12. PACKING LIST

Principal machine:1
 Nasal cannula:2
 Humidifier connection tube:2
 Primary air suction filter:3
 Secondary air suction filter: 8
 Plastic adapter for secondary filter: 2
 User Manual: 1

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Fault symptom	Check	Measures
After switching on or during operation, the LCD display does not show anything and an alarm sound is heard.	1. Power failure? 2. Is the power cable loose?	1. Use the equipment after switch-on. 2. Check and insert the power cable correctly.
Purity indicator turns yellow or red.	1. Is the oxygen flow too large? 2. Does the molecular sieve age after prolonged use?	1. Adjust the oxygen flow. 2. Replace the molecular sieve.
	Note: Replacement of the molecular sieve should be carried out by professionals	
Oxygen concentrator operates normally, but no oxygen is coming out of the tubing.	1. Has the oxygen tube been inserted correctly? 2. Is the humidifier bottle cap tight?	1. Adjust the oxygen flow. 2. Replace the molecular sieve.
Oxygen concentrator operates normally, but no oxygen is coming out of the tubing.	1. Has the oxygen tube been inserted correctly? 2. Is the humidifier bottle cap tight?	1. Insert the oxygen tube correctly. 2. Tighten the cap of the humidifier bottle.
The equipment works normally, but the flow cannot be adjusted to the maximum value, the bubbles in the humidifier bottle are small.	1. Is the air suction filter clogged with dust? 2. Is the oxygen tube bent?	1. Clean or replace the filter. 2. Unbend the oxygen tube.
Noisy during operation.	Has the concentrator been installed without any problems?	Install the equipment according to the recommendations.
Fault alarm indicator turns red.	Low voltage.	It is recommended to install a voltage stabiliser.
Fault alarm indicator turns red.	Circulation pressure is low.	1. Clean or replace the filter. 2. If this situation continues, please contact your distributor.
Fault alarm indicator turns red.	Circulation pressure is high.	Stop using the device immediately, contact your dealer.
Fault alarm indicator turns red.	Air compressor failure.	Stop using the device immediately, contact your dealer.
No alarm sound during a blackout and the alarm indicator does not light up.	Is the battery dead?	Replace the battery.
Other faults.	Contact the distributor.	

Check if the alarm system works

To check if the alarm system is working, disconnect the power cord when the oxygen concentrator is operating normally. If the fault indicator turns red and the alarm sound is emitted, then the alarm system is operating normally, otherwise contact the distributor for inspection or repair.

8. MAINTENANCE

8.1 Maintenance considerations

■ Only dealers authorised by the manufacturer or qualified personnel after special training may perform preventive maintenance or performance adjustments.

In case of malfunctions contact the distributor and ask for inspection/repair of the equipment; do not disassemble the equipment for repair yourself.

■ Clean the humidifier bottle every day. Add pure (distilled) water to the humidifier bottle and the water level should be between the maximum and minimum levels.

Regularly clean the air suction filter foam and filter cotton. Clean beforehand in case of excessive dust.

■ The oxygen concentrator must not be started before installing the filter cotton and filter foam.

■ The power cable is supplied according to the power of the equipment; do not replace it arbitrarily.

8.2 Cleaning of the humidifier bottle

■ To clean, first remove the silicone rubber tubing from the humidifier bottle, then remove the humidifier bottle. The humidifier bottle is usually cleaned with water. Pay attention to clean the small air inlet hole at the bottom of the central tube in the bottle, to keep the oxygen flow unobstructed.

■ After cleaning, sterilise (soak in disinfectant containing chlorine at a concentration of 500 mg/L) and rinse with distilled water. Add the appropriate amount of pure (distilled) water to the bottle; the water level should be between the maximum (MAX) and minimum (MIN) lines; then tighten the bottle cap and place it on the transparent cover and insert the silicone rubber tubing.

8.3 Cleaning of primary air suction filter (filter sponge) and secondary filter

■ Timely cleaning and replacement of the primary air suction filter and secondary filter is very important to protect the air compressor and molecular sieve, and to extend the life of the oxygen concentrator.

■ The primary air suction filter should be cleaned once a week. Remove the filter sponge and rinse it with clean water. It can be reused after air-drying.

■ The secondary air suction filter shall be cleaned once every half a month. Unscrew the filter clockwise, take it out and rinse it with clean water. After it has air-dried, install it in the original place and tighten it counterclockwise.

8.4 Cleaning the oxygen concentrator

■ After switching off the equipment, use a soft towel moistened with a small amount of neutral household cleaner to wipe the housing, and then use a dry towel to dry. Be careful to prevent liquid from seeping into the slot in the housing.

8.5 Battery replacement

When the sound of the power failure alarm is significantly reduced (user can verify this by disconnecting the power supply), the battery needs to be replaced. Battery specifications: 9V alkaline battery. To replace, open the battery compartment cover on the transparent cover, remove the exhausted battery, insert the new battery according to the marked polarity.

9. TROUBLESHOOTING

In case of problems during use, before sending the equipment for repair, read this section carefully; you may be able to resolve the situation yourself. If the problems cannot be solved according to the examples in Table 2, send the oxygen concentrator for repair; do not attempt to repair it yourself.

INDEX OF LANGUAGES

Spanish	2-12
English	13-22
French	23-30

INDEX OF CONTENTS

1. Brief description of the product	13	7. Product use	17
2. Safety overview.....	13	8. Maintenance	20
3. Principle of operation.....	14	9. Troubleshooting	20
4. Structural features.....	15	10. Transportation and storage	22
5. Technical parameters.....	15	11. Accessories	22
6. Product installation	16	12. Packing list.....	22

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The AXAVET OLV-5 oxygen concentrator uses air as raw material and high quality molecular sieve as adsorbent, it adopts the principle of pressure swing adsorption (PSA) to directly separate oxygen from nitrogen at normal temperature, thus producing high purity oxygen.

To ensure the safety and efficiency of the equipment read this manual carefully before use, in order to have a thorough understanding and knowledge of the product's performance, as well as the correct methods of operation and maintenance. Strictly observe the relevant safety precautions during installation, use and maintenance.

1.1 Product features

■ Plastic housing, novel design, simple and stable operation, easy to maintain.

■ It generates oxygen at low cost by adopting physical methods, with air as raw material, without the use of additives, only needing an electrical power supply.

Adopts efficient molecular sieve pressure swing adsorption (PSA) technology, with simple process flow and low energy consumption.

2. SAFETY OVERVIEW

2.1 Safety precautions

■ This product cannot be used as life support for critically ill patients who require additional care.

Oxygen therapy may be harmful under certain conditions. The health care professional must monitor oxygen flow and absorption time.

■ Patients with severe carbon monoxide poisoning should not use this product.

■ In anticipation of possible power outages or malfunction of the oxygen concentrator, other backup oxygen supply devices (e.g. oxygen cylinder, oxygen bag, etc.) should be available.

- If the oxygen concentrator is not functioning properly, discontinue use immediately and seek technical assistance from the supplier to resolve the problem.
- Oxygen is a gas involved in combustion, therefore this equipment must not be used where open flames or flammable substances are present.
- The power supply must comply with electrical safety standards. The oxygen concentrator shall not be used when the protective earth terminal of the power supply does not comply with the relevant regulations, otherwise personal injury may occur.
- Turn off the power and unplug the power cord before cleaning and maintaining the oxygen concentrator.
- Unauthorised persons shall not open the equipment for maintenance.

2.2 Guidance on the electromagnetic environment


- This product is suitable for use in buildings connected to the civil low-voltage mains power supply.
- The RF energy used by this product is for internal operation only. Therefore, its RF emission is very low, with no impact on other electrical equipment in the vicinity.
- Due to the impact of radio transmitting devices or other sources of electrical noise in health care facilities, severe interference caused by too close a distance or high transmitting power may cause this product to be disrupted.
- If this happens, check the places of use to find out the source of interference and take the following measures to eliminate the interference: 1) Turn off nearby devices before turning on the machine; 2) Change the direction or location of interfering devices; 3) Increase the distance between interfering devices and this product.

2.3 Environmental protection

Waste disposal must comply with local legal regulations.

3. PRINCIPLE OF OPERATION

The process flow of the oxygen concentration process is shown below:

Alarm	Cause	Audible alarm	Color of indicators	LCD display
Oxygen concentration alarm	$65\% (\pm 3\%) \leq \text{Oxygen concentration} < 82\% (\pm 3\%)$	Yes	Purity indicator in yellow	N/A
	Oxygen concentration failure $\leq 65\% (\pm 3\%)$	Yes	Purity indicator in red	Low concentration
High pressure alarm	Circulating system pressure is high	Yes	Trouble indicator in red	Shutdown High pressure
Low pressure alarm	Circulating system pressure is low	Yes	Trouble indicator in red	Low pressure
Air compressor failure alarm	Short circuit of compressor	Yes	Trouble indicator in red	Shutdown Air compressor failure

■ Power failure alarm

When the mains voltage is suddenly interrupted (power failure), the equipment will stop working, emit an alarm sound and the failure indicator will turn red.

■ Low voltage alarm

When the mains voltage falls below $180V \pm 10V$, the oxygen concentrator will stop working, sound an alarm, the fault indicator will turn red and the LCD display will indicate "Low voltage". It is recommended to install a voltage stabiliser to prevent the occurrence of this fault.

■ Oxygen concentration alarm

After the concentrator starts and generates oxygen for 10 minutes, under the circumstance that the output oxygen flow is not greater than the maximum recommended flow:

When the oxygen concentration is $82\% (\pm 3\%)$ or more, the oxygen purity indicator is green. The equipment is operating correctly.

When the oxygen concentration is within the range $65\% (\pm 3\%) < 82\% (\pm 3\%)$, the oxygen purity indicator is yellow and an alarm sound is emitted. Contact the distributor immediately. The user may temporarily continue to use the machine if another source of backup oxygen is available.

When the oxygen concentration is below $65\% (\pm 3\%)$, the purity indicator is red, an alarm sound is emitted and the LCD display indicates "Low Concentration". Immediately turn off the equipment. Use the backup oxygen source. Contact the distributor for inspection or repair.

■ High/low pressure alarm

When the circulating system pressure is high, the equipment shuts down, emits an alarm sound, the fault indicator is red and the LCD display indicates "High pressure".

When the circulating system pressure is low, the concentrator emits an alarm sound, the fault indicator is red and the LCD display indicates "Low pressure".

In case of high or low pressure alarm, use backup oxygen and contact the distributor for inspection or repair.

■ Air compressor failure alarm

In case of short-circuit of the air compressor, the unit shuts down, gives an alarm sound, the fault indicator is red and the LCD display shows "Air compressor fault". Use reserve oxygen and contact the dealer for inspection or repair.

■ **Add water to the humidifier bottle:** Remove the silicone rubber tube from the bottle, take out the bottle. Add the appropriate amount of pure water to the bottle (distilled water), the water level should be between the MAX and MIN lines. Tighten the cap of the bottle, place it in position and insert the silicone rubber tube.

■ **Switching on:** Connect the power cable, press the ON/OFF switch, the “ON” indicator will light up. The equipment is ready to start operating.

■ **Continuous / timed oxygen generation:** Press the “oxygen generation / timing” button, the concentrator will turn on and enter the “continuous oxygen generation” state. The LCD display shows the total usage time and other information, as shown in Fig. 5. Press the “oxygen generation / timing” button again, the concentrator will enter the “timed oxygen generation” state and the time is 30 minutes (see Fig. 6); for each press of the “oxygen generation / timing” button, the preset time will increase by 30 minutes (the maximum preset time is 120 minutes); the user can preset the time according to his requirement. If after setting a time of 120 minutes the “oxygen generation / timing” button is pressed again, the device will stop generating oxygen.

Total time: 00000 H
Oxygen is being generated
30 minutes

Fig. 6 Display of timed oxygen generation

■ **Adjust the flow rate:** Open the transparent cover, slowly turn the flow control knob in the marked direction until the desired value is reached (the reading will be referenced to the centre of the black float). Meanwhile, there are bubbles in the water in the humidifier bottle and there will be oxygen supply from the “oxygen outlet”.

■ **Start oxygen absorption:** Insert the oxygen tubing into the “oxygen outlet” of the humidifier bottle, close the transparent cover. After confirming that there is oxygen output, use the tubing (refer to the oxygen tubing operating instructions for methods of use) to start oxygen absorption.

■ **Automatic shutdown:** After reaching the preset time, the equipment will shut down automatically. The LCD display is shown in Fig. 7.

Total time: 00038 H
Shutdown
Time out

Fig. 7 Shutdown display

■ **Immediate shutdown:** During operation, press the ON/OFF switch on the top cover to immediately turn off the equipment.

7.3 Alarm functions

After the concentrator starts and generates oxygen for 10 minutes, during normal operation, the “Power”, “Purity” and “Trouble” indicators will be green.

If the faults listed in Table 1 occur, the equipment will activate an audible / visual alarm. The alarm functions are shown in Table 1.

Alarm	Cause	Audible alarm	Color of indicators	LCD display
Power failure alarm	Power supply interruption	Yes	Trouble indicator in red	N/A
Low voltage alarm	Power supply voltage is less than 180V±10V	Yes	Trouble indicator in red	Low voltage

The oxygen concentrator uses molecular sieve as adsorbent, adopts the principle of pressure swing adsorption (PSA), sends the air after filtration to the molecular sieve adsorption column for the cyclic process of pressurised adsorption and decompressed desorption; thus high purity oxygen is produced.

4. STRUCTURAL FEATURES

5. TECHNICAL PARAMETERS

5.1 Environmental conditions

Ambient temperature: 10 °C ~ 40 °C

Relative humidity: 30% ~ 75%

Atmospheric pressure: 860hPa ~ 1060hPa

No corrosive gases or strong magnetic fields must be present in the vicinity of the equipment.

5.2 Inlet air requirements

- Impurities: ≤ 0,3 mg/cm³

- Oil content: ≤ 0,01ppm

5.3 Product functions

Time function: shows the total working time on the display screen

Time setting: the required oxygen absorption time is set

- Automatic shutdown: the equipment shuts down automatically after reaching the preset oxygen generation time
- Power failure alarm function
- Low voltage alarm function
- Low oxygen concentration alarm function
- Circulation pressure failure alarm function
- Air compressor short-circuit alarm function

5.4 Technical indicators

- Physical and chemical indicators of oxygen
 - Oxygen concentration: $\geq 90\%$ (V / V)
 - Water content: $\leq 0,07 \text{ g/m}^3$
 - Carbon dioxide content: $\leq 0,01\%$ (V / V)
 - Carbon monoxide content: meets the requirements of table 1 in GB 8982-2009
 - Acid and alkali gas content: meets the requirements of table 1 in GB 8982-2009
 - Ozone and other gaseous oxides content: meets the requirements of table 1 in GB 8982-2009
 - Odour: odourless
 - Particle size of solids: $\leq 10 \mu\text{m}$
 - Solids content: $\leq 0,5 \text{ mg/m}^3$
- Technical specifications of product
 - Recommended maximum flow rate (when oxygen concentration $\geq 90\%$): 5 L/min
 - Flow range when the nominal outlet pressure is 7 kPa (oxygen concentration $\geq 90\%$): 0-5 L/min
 - At maximum recommended flow rate, flow rate changes when 7 kPa back pressure is applied: $< 0.5 \text{ L/min}$
 - Flow range when nominal outlet pressure is 0 (oxygen concentration $\geq 90\%$): 0-5 L/min
 - Noise level during operation: $< 45 \text{ dB}$
 - Timer error: no more than 3%
 - Output pressure: 100-115 kPa
 - Air compressor safety valve release pressure: 200-300kPa
 - Oxygen concentration measurement: 0-95%, error $\leq \pm 3\%$
 - Mode of operation: continuous/timed operation
 - Input power: 500VA + 15%
 - Power supply: AC 220V $\pm 22\text{V}$, 50Hz $\pm 1\text{Hz}$
 - Weight: 22 kg
 - Dimensions (L x W x H): 490x270x560mm

6. PRODUCT INSTALLATION

6.1 Unpacking inspection

Open the cardboard box from the top; then open the plastic bag and lift the oxygen concentrator by grasping the front and rear handles. Carefully check for any damage during transport and also the accessories and documents according to the packing list.

6.2 Installation precautions

- The oxygen concentrator shall be installed in ventilated indoor locations, free of dust, corrosive, toxic or harmful gases or fumes. Avoid direct sunlight. The distance from walls and other objects should be greater than 10 cm.
- The oxygen concentrator must not be installed where there is an open flame, source of fire, danger of explosion, high humidity, too high or too low temperatures. Furthermore, it must not be used in an enclosed space.
- Sundries, water or oil containers shall not be placed on top of the equipment.

- The concentrator must not be placed on soft surfaces (e.g. beds, sofas) that may cause tilting or sagging. Avoid shutdown or decrease in oxygen concentration caused by too high a temperature due to blockage of the air inlet or outlet.
- The oxygen concentrator must be positioned carefully, otherwise it will increase the noise during operation.
- If the mains voltage is unstable and outside the range of 220 $\pm$ 22V, install a voltage stabiliser before using the equipment.
- Be sure to install the battery before use, otherwise some alarm functions will be lost.
- Connect the concentrator to a safe and qualified socket.

7. PRODUCT USE

7.1 Use precautions

- During use, make sure that the exhaust is not obstructed at the bottom of the concentrator, otherwise internal overheating may occur.
- When the oxygen output is less than the maximum recommended flow, the oxygen concentration reaches 90%. When the flow exceeds the maximum recommended flow, the oxygen concentration will decrease with increasing flow.
- The concentrator will reach the specified performance after 10 minutes of operation.
- There will be an intermittent exhaust sound (every 6 seconds) during operation, which is normal.
- No oil, grease or other similar substances shall be used on or near the concentrator, and no lubricants other than those recommended by the manufacturer shall be used.
- During use, add water in a timely manner when the water level in the humidifier bottle is below the minimum level.
- The concentrator must not start immediately after a shutdown; wait 5 minutes to restart it.
- When the indicated oxygen concentration is abnormal, stop using the equipment and contact the distributor for inspection and maintenance.
- The molecular sieve will age due to use, environment, etc., causing an irrecoverable decrease in the amount of oxygen generated. When this phenomenon occurs, contact the distributor to replace the molecular sieve.
- If the concentrator will be idle for a long period of time, be sure to unplug the power cord.

7.2 Operation

- The control panel of the concentrator is shown in Fig. 3.