M A N U A L D E U S O

Microscopio profesional

Serie 158

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Operation Manual
For
Research Biological Microscope
Series 158

In order to exert performance of this microscope completely, please read Operation Manual carefully before use and make it well-kept after use.

Directory

Prologue	1
I Specification	2
II Configuration	3
III Installation	7
IV Detailed observation steps	8
1. Brightness adjustment	8
2. Specimen emplacing	8
3. Focusing	8
4. Drawtube adjustment	9
5. Interpupillary distance adjustment	9
6. Diopter adjustment	10
7. Adjusting condenser position & aperture diaphragm	10
8. Mounting filters	10
9. Field diaphragm adjustment	11
10. Changing objectives	11
11. Using 100X oil immersion objective	11
12. Light source adjustment	12
13. Changing the bulb	12
14. Connecting power cable	13
15. Replace the fuses	13
V Problems occurred in using and disposal method	14

Prologue

Series 158

Research Biological Microscope

The diversification of this model brings you enter into an epoch not only the performance, but also the application. The module design can meet your needs on all normal experiment and scientific research. Its humanized configuration design, reasonable control device and layout of every operation knob can keep your eyes, neck, shoulders, waist, arms and wrist against tiredness. You will be released from heavy work due to the comfortable operation.

Important advice

- 1) Don't let the microscope exposure under Sun light directly, don't placed it in where is high temperature、wet、dust or acid and alkali.
- 2) The microscope must be earthed in order to avoid getting an electric shock.
- 3) The power switch must be closed and cut off the electrical power when changing the bulb or fuse. Please don't change the bulb until the lamp is cool completely.
- 4) You must use specified bulbs with the same specification supplied by our company.
- 5) Don't knock down any part of the microscope in avoid of damage and trouble.
- 6) After finishing work, please cut off the power supply and cover the dustcover.
- 7) Clean the lens with gauze moistened with a spot of mixture of ether and alcohol 7:3.

Safety Symbol

Mark	Description
	Hot surface, don't touch Scalded danger
	Read Operation Manual before use. Improper operating will cause user injured or microscope damaged
I	Main power source Turn On
0	Main power source Turn Off

I Specification

Optical system	Corrected Infinite optical system
Viewing head	Compensation free binocular head, Incline adjustment $0^{\circ} \sim 40^{\circ}$ Compensation free trinocular head, Incline adjustment $0^{\circ} \sim 40^{\circ}$ Beamsplitter rate: Binocular head 50%, Trinocular head 50%
Interpupillary distance	55mm~75mm
Eyepiece	WF10×/22mm × ()
Nosepiece	Quintuple nosepiece ()
Objective	Infinite Plan Objective: 4x, 10x, 20x (S), 40x (S), 100x (S)Oil×, ×, ×× ×, × ()
Mechanical stage	Stage: 180mm×160mm Moving range: 75mm×50mm Wide mechanical stage with 2-Microslide specimen clamp
Condenser	Swing type N.A.1.25 Abbe condenser with iris diaphragm and filter holder
Focusing	Coaxial coarse & fine focusing adjustment with rack and pinion mechanism Fine focusing scale value 0.002mm
Light collector	Kohler illumination
Illumination	Halogen lamp 12V/50W Adjustable brightness
Optional accessory	fluorescent kit Phase Contrast unit Dry or Immersion dark field condenser Simple polarizing kit 3.0 or 1.3 Mega pixels CMOS electronic eyepiece
Operating environment	Indoor use Temperature: 5℃~40℃ Maximum relative humidity: 80%

Objectives

Feature Objective	Magnification	Numerical aperture (N.A.)	Working distance (mm)	Cover glass thickness (mm)	Remark
Infinite plan objective	4×	0.10	29.40	0.17	
	10×	0.25	16.01	0.17	
	20×(S)	0.40	3.00	0.17	
	40×(S)	0.65	0.63	0.17	
	60×(S)	0.85	0.50	0.17	Optional
	100×(s)Oil	1.25	0.35	0.17	

II Configuration

—158/7


Fig. 1.

—158/7


Fig 2.


Fig. 3


Fig. 4

III Installation


Fig. 5.

1. Objective
5. Power source cable

2. Swing type condenser
6. Humanized binocular head

3. Base
7. Eyepiece

4. Lamp housing

IV Detailed Observation steps

1. Brightness adjustment (Fig. 6.)


Fig. 6.

- 1) Press main power switch ① to “I” side.
- 2) Revolving brightness knob ② as per arrowhead to brighten and darkle rightabout. The arrowhead on the switch panel shows voltage value.

2. Specimen emplacing (Fig. 7.)


Fig. 7.

- 1) Revolve coarse focusing knob ② to lower the Mechanical stage according to arrowhead direction
- 2) Open the spring clamp ③ as per arrowhead direction and put the specimen on the stage from the front.
- 3) Clip the specimen by laying back spring clamp ③
- 4) Revolve the knob ④ to make specimen longitudinal movement and knob ⑤ for transversal movement.

★ Please don't push mechanical stage by hand directly to move the specimen in avoid damage of transmission parts.

3. Focusing (Fig. 8.)


Fig. 8.

Main points of focusing

- 1) First try to make the object close to specimen, then use 10X objective and rotating coarse focusing knob ① as per arrowhead.
- 2) After finding the specimen, please make fine focusing by adjusting knob ②


Fig. 9.

Adjusting tension of coarse focusing knob (Fig. 9.)

- 1) Tension of coarse focusing knob is adjustable. Adjusting the ring ① as per arrowhead direction to make heavy tension and reverse way to make light tension
- 2) If you find the mechanical stage slides down automatically or the image fogs immediately after fine focusing, it is caused by loose of tension ring, please fasten the tension ring as per arrowhead direction

4. Drawtube adjustment (Fig.10.)

Fig. 10.

Please adjust an angle from $0^{\circ} \sim 40^{\circ}$ of the drawtube to meet your requirement, and try to make you comfortable under long time working.

5. Interpupillary distance adjustment (Fig. 11.)

Fig.11.

◎ Adjusting the interpupillary distance between your 2 eyes, and make them aim to one image under the eyepiece at one time. That will prevent tiredness during observing. Method: observing the image through 2 eyepieces and adjusting 2 drawtubes to make viewing field consistent both left and right.

6. Diopter adjustment (Fig. 12.)


Fig. 12.

◎ Adjust eyesight of left & right eye of observer

- 1) Firstly observe right side eyepiece by right eye, rotating coarse & fine focusing knob till you see clear image of specimen
- 2) Then observe left side eyepiece by left eye, rotating Diopter ring ① to also make specimen image clear

Fig. 13.

Usage of eyeshade (Fig. 13)

For people wear glasses

The eyeshade can be folded (prevent glasses touch eyepiece directly)

For people without glasses

Eyeshade can prevent interreflected incident light between eyes and eyepiece, making observation easily.

7. Adjusting condenser position & aperture diaphragm (Fig. 14.)


Fig. 14.

◎ In general, the condenser is at upper stop position, but when observing the field brightness which is not uniformity, please adjust the condenser a little down to get better illumination.

- 1) Moving condenser to the top by adjusting knob ①
- 2) Adjusting rod of iris diaphragm of condenser to make the aperture diaphragm consistent with numeral aperture of objective to get clear and sharp image. adjusting rod ② of iris diaphragm of condenser.
- 3) When use 4X objective, please insert the swing type condenser as per arrowhead direction to fill the field with full light.

8. Mounting filters (Fig. 15.)

Fig. 15.

◎ Swing out the filter holder, put the filter in it and then backtrack

9. Field diaphragm adjustment (Fig. 16.)

Fig. 16.

Adjusting Field diaphragm ring to make illumination area fitted with eyepiece field.

If the Field diaphragm is bigger than aperture required, external light will come into the viewing field, that will interfere the quality of image and the contrast

10. Changing objectives (Fig. 17.)

Fig. 17.

Rotating nosepiece ① to make required objective over the specimen, make sure the nosepiece is clicked on correct position.

11. Using 100X oil immersion objective (Fig. 18.)


Fig. 18.

◎ When use 100X (oil) objective, please drop oil on the inswept part of objective, otherwise, clear image can not be obtained.

- 1) Focusing sequence is from low magnification objective to high magnification
- 2) Please drop oil ① on the observing point of the specimen before 100X (oil) objective shift to the light path
- 3) Rotating nosepiece to turn 100X(oil) objective into light path, then adjust fine focusing knob to obtain clear image

★ If there's air bulb in oil, it will influence observation. Take out air bulb by swinging nosepiece several times

- 4) Wipe up the oil on objective and specimen carefully by gauze moistened with a spot of 7:3 ether and alcohol mixture

12. Light source adjustment (Fig. 19.)


Fig. 19.

Adjusting position of filament to obtain better illumination

- 1) Filament forward or backward adjusting knob ①
- 2) Filament up or down adjusting knob ②
- 3) Filament left or right adjusting knob ③

13. Changing the bulb (Fig. 20.)


Fig. 20.

★ The electrical power must be disconnected when changing the bulb.

- 1) Unscrew the screw ① on the top of the left side, hold the cover and take it off softly
 - 2) The halogen bulb ② should not have fingerprint on it, don't touch bulb by finger directly when changing the bulb, please wear gloves or encase the lamp with clean paper when inserting bulb into socket.
- ◎ The fingerprint or dirty on the bulb will reduce the brightness, shorten the life of the bulb, please clean it with alcohol moistened cloth
- 3) Close the lamp cover after the bulb is changed

14. Connecting power cable (Fig. 21. Fig. 22. & Fig. 23.)


Fig. 21.

- 1) Turn the power switch ① to “O” (OFF), then connect the power cable.
- 2) Insert the plug ② into the inlet ③, and make sure the connection is good.
- 3) Insert another end plug ④ of the power wire to the power supply socket ⑤

15. Replace the fuse (fig.22)


Fig. 22.

- 1) Turn the main power switch ① to “O” (OFF) before replacing the fuse.
- 2) Take out the fuse ⑥ from the fuse box ⑦, then replace it by new one


Fig. 23.

V Problems occurred in using and disposal method

According to different use method, although it is not malfunction, the microscope performance may not be given out completely. If the trouble occurred, please try to find the cause.

Problem point	Cause	Disposal method	page
1. Brightness of viewing field is not uniformity	Objective is not entered into the light path	Ensure the objective is entered into the light path ("click" sounded)	11
	Condenser position is too low	Adjust to the top	10
	Dirty or dust on objective, eyepiece, condenser and light exit	Cleaning completely	1
2. You can see the dirty and dust in the viewing field	There's dirty and dust on eyepiece condenser, light exit and the specimen	Cleaning completely	1
3. Viewing image glaring	Condenser is too low	Drive up the condenser	10
4. Viewing image looks frosted and unclear	Objective is not entered into the light path	Ensure the objective is entered into the light path ("click" sounded)	11
	There's dirty and dust on eyepiece, objective, condenser, light exit and the specimen	Cleaning completely	1
	there's no immersion oil on oil immersion objective	Using immersion oil	11
	Air bulb in immersion oil	Take out the air bulb	11
5. When focusing with high magnification objective, it will touch the specimen	The specimen placed is on the reversed side	Reversing Cover glass side on the top, then place the specimen again	8
6. Focusing adjusting knob is too heavy	Tension adjusting ring is too tight.	Loosen properly	9
7. The stage goes down automatically, the image becomes blurry during observing	Tension adjusting ring of the coarse focusing adjusting knob is too loose.	Tighten properly	9
8. The coarse focusing adjusting doesn't go down	Condenser is too low	drive up the condenser	10
9. Field of vision of 2 eyes differed	interpupillary distance is not fitting	adjusting correctly	9
	Diopter of 2 eyes is not adjusted equally	adjusting correctly	10
	Magnification of left and right eyepiece is different	Change to the same magnification eyepiece	2
10. Objective touches specimen, while change the objective from the low magnification to high magnification	The specimen placed is on the reversed side	Reversing Cover glass side on the top, then place the specimen again	8
	The cover glass is too thick	Change to cover glass 0.17mm thick	8
11. The bulb does not work	Forget to plug in	Make sure to plug in	13
	The bulb is not inserted correctly	Insert it correctly	12
	The bulb is burn out	Replace bulb with a new one	12
12. The field of view is not bright enough	Brightness adjustment knob turned to wrong position	Adjust correctly	8
13. The bulb light flickers or the brightness is not stable	The bulb uses long time and will be burn out soon	Replace it by a new one	12
	The wire is not connected in order	Connecting wire completely	13

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© Ajuste la distancia interpupilar entre los 2 ojos de modo que consiga ver una única imagen con ambos ojos a la vez. Esto evitará fatigas durante la observación.