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APPARATUS FOR HEAT CONDUCTIVITY OF METALS

Article code: QLJ002

Parts Included:

1. Conductivity of metals apparatus
2. Small quantity of wax
3. Instructions

Additional Items Required:

1. Heating device such as alcohol lamp or Bunsen burner
2. Laboratory stand
3. Timing device such as stop-watch

Theory:

When any substance is heated in a small portion, temperature of that portion rises. This rise in temperature is gradually transmitted to other portions of the substance till temperature equilibrium is established. This can be observed by heating a rod of a material at its one end while holding it in the hands from the other end. After some time, a sensation of hotness is felt. Thus, it is evident that, in an object, heat travels from one part of the object at higher temperature to its other parts at lower temperature. The phenomenon of transmitting heat in this manner is possessed by all the substances to a varying degree and is called conduction of heat. Accordingly, the substances are classified into good conductors of heat and bad conductors of heat or insulators. In case of conductors, such as metals, there is a large number of free electrons actually separated from the atoms and hence capable of moving freely in the inter-atomic spaces. The thermal conduction phenomenon in metals is assumed to be affected solely by the motion of free electrons which transmit thermal energy from one end having higher energy (or higher temperature) to that of lower energy (or lower temperature) in the form of transmission of kinetic energy. Electrons at the hot end possess more kinetic energy and this kinetic energy of electrons decreases gradually on moving away from hot end. These electrons, having higher energy, transfers kinetic energy to the others having lower energy through collision, thus progressively increasing the energy of electrons away from hot end, resulting in the increase of their temperature. In this way, heat is conducted from hot end to the cold end of metallic conductors by the process of conduction.

As discussed above, the rate of conduction of heat from one end to other end of a metallic rod with respect to time (also called thermal conductivity of the metal) depends upon the availability of number of free electrons per unit

volume. Since different metals have different number of free electrons available per unit volume, their conductivity also varies.

In case of insulators or poor conductors of heat, the electrons are so strongly bound to the atom that they can not move so freely as compared to that of conductors except for a very small strained displacement about the mean equilibrium position, which is usually oscillatory in nature.

Description of the Apparatus:

The apparatus consists of a hexagonal metal block fixed at one end of a metallic rod. An insulated handle is fitted at the other end of the rod. Five metal rods of same length and diameter, but each of different material are mounted radially on the hexagonal metal block (as shown in the figure). These five rods consist of copper, aluminum, brass, iron and stainless steel. Each of the five rods has a small cavity at its free end to hold small lump of wax. Relative conductivity of these five metal rods can be studied.


Figure: Conductivity of Metals Apparatus

Experiment:

1. Put a small lump of wax in the cavity at the free ends of each of the five metal rods.
2. Hold the conductivity apparatus in a clamp on a laboratory stand. If stand is not available, it can also be safely held from the insulated handle in your hands.

3. Light-up an alcohol lamp or Bunsen burner and position the flame below the central hexagonal metal block.
4. The heat is conducted through each of the metallic rods. By observing the sequence in which the wax put in each of the rods melts, their conductivity can be compared. The rod on which wax melts first is the most conducting and on which the wax melts last is the least conducting of all the metals.
5. Time taken by each rod to completely melt the wax can be observed using a timing device such as stop-watch, which gives an indication of the degree of variance in their conductivities.

It will be observed that copper is most conducting and stainless steel is least conducting of all the five metals used.

Precautions:

1. While using alcohol lamp or Bunsen burner proper safety regulations and cautions must be observed to avoid any accident.
2. The apparatus should be held firmly from the handle with all the five rods positioned along horizontal plane.
3. Ensure the same level of flame for uniform heating of the metal rods with respect to time for better experimental results.
4. Let the apparatus cool down properly before touching any of the metal parts.