Imán de herradura, acero cromado

Referencia QLL004

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Horseshoe magnet, chromed steel

Code QLL004

Horseshoe magnets are just bar magnets bent in a U shape. The U shape makes the magnet stronger by pointing the poles in the same direction. Originally created as a replacement for the bar magnet, this shape has become the universal symbol for magnets. It can be used to pick up metal objects of any size depending on the strength of the horseshoe magnet. For example, small horseshoes can collect paper clips while industrial-sized horseshoe magnets are used in construction and engineering to pick up large pieces of heavy metal. Horseshoe magnets are also used at the bottom of pendulums.

Experiment 1: To observe the magnetic field lines with the help of iron filing and compasses.

Components Required:

- Horseshoe magnet
- Compasses (not included)
- Large sheet of white paper (not included)
- Iron filing (not included)

Procedure:

1. Place a large white sheet on the table.
2. Now place a horseshoe magnet on the sheet.
3. Sprinkle iron filing on the magnet, the iron filing get attracted towards the ends of the magnet and form the pattern according to the magnetic field lines of the magnet.
4. Place the compasses near the magnet to observe the direction of the magnetic field lines.

Horseshoe magnet