DEMOSTRADOR DE CAMPO MAGNETICO

MAGNETIC FIELD DEMONSTRATION

REF: QLL008 -RECTO – STRAIGHT

REF: QLL009- BUCLE CIRCULAR COIL

REF: QLL10 -BOBINA SOLENOID

INDEX OF LANGUAGES- ÍNDICE DE IDIOMAS

ESPAÑOL	2
ENGLISH.....	3

DEMOSTRADOR DE CAMPO MAGNETICO

Descripción

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
MAGNETIC FIELD DEMONSTRATION

Description

For the study and demonstration of magnetic fields associated to different shapes of an electrical current conducting wire. Mounted on a transparent acrylic base they can be placed on a retro projector for classroom demonstrations. Each instrument presents a pair of sockets for 4 mm banana plugs to be connected to a power supply. The effects of the magnetic field can be observed by using iron filings.

3 models are available depending on the shape of the conducting wire

1. Straight: U-shaped copper wire vertically mounted on the base. Maximum current of 8 A.
2. Circular coil: 1 coil copper wire. Maximum current of 8A.
3. Solenoid: copper solenoid of about 52x125 mm diameter x length

Operation

1. Connect the power supply with the solenoid through the banana sockets
2. Scatter iron filings on the stand of the solenoid in such a way that only the approximate surface around the solenoid is covered.
3. Switch on the power supply and gently tap the stand
4. You will observe the ring shaped magnetic lines of force around the solenoid.
5. As soon as you observe the magnetic lines of force, switch off the power supply.

In this experiment using solenoid, several current carrying conductor loops are added together. According to Fleming's right hand rule if the thumb of your right hand is pointing up in the direction of the current then your fingers point in the direction of the magnetic field.