Limadura de hierro

Referencia QLL005

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Iron filings

Code QLL005

Iron filings are very small pieces of iron that look like a light powder. They are very often used in science demonstrations to show the direction of a magnetic field. Since iron is a ferromagnetic material, a magnetic field induces each particle to become a tiny bar magnet.

Here are some experiments where we can use iron filings.


Experiment 1: To observe the magnetic field lines of a bar magnet and observe the magnetic field lines using compasses.

Other Components Required:

- Bar Magnet
- Six compasses
- Large white paper

Procedure:

1. Place a large white sheet on the table.
2. Now place a bar magnet on the sheet.
3. Sprinkle iron filling on the bar magnet; the iron filing get attracted towards the end of the bar and form the pattern according to the magnetic field lines of the bar magnet.
4. Place the compasses near the magnet to observe the direction of the magnetic field lines.


Experiment 2 : To determine the pole of the natural bar magnet.

Components Required:


- Natural bar magnet pair
- Large white paper
- Round shape magnet pair

Procedure:

1. Place a large white sheet on the table.
2. Now place two natural bar magnet such that they repel each other.
3. Now sprinkle some iron filing over the magnets .
4. Iron filings spread according to the magnetic field of both the magnets.

5. Again place the magnets so that they attract each other and repeat the same procedure.

Again repeat the same experiment using round magnets.