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Autoclave's qualified installation

On the EC area, Directive 2014/68/UE provides for the qualified installation of this autoclave.

This installation should be performed by a "qualified installer of pressure vessels".

You may request information from your dealer for a list of qualified installers in your geographical area.

Documents needed to handle the qualified installation are the following:

- EC Declaration (document that comes with the device).

Some other documents may also be obtained from **service@jpselecta.es**:

- Test run.
- EC Declaration of Safety Valve.
- The qualified installer also will provide you information about the periodic inspections.

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1 Information for users

1.1 Autoclave incoming

Handle the equipment with care. Unpack and check that all items coincide with the delivery note. If you see any discrepancy or damage of any kind, notify your distributor.

Keep original package for some days. After this time, dispose of it properly, separating its components: cardboard, wood and plastic.

1.2 Documentation

This documentation comes with the package:

- Instruction manual

Read this manual before using the equipment.

These instructions are an important part of the equipment and should be made available to all users.

If you are in any doubt or want to clarify the use of this instrument, contact your local or the service department of J.P. Selecta, s.a.u.

- EC declaration

1.3 Other important information

The autoclaves sent to J.P. SELECTA, s.a.u. for repairs or maintenance must be decontaminated. In case of failure to comply with this rule, the equipment will be returned to the sender.

Any unauthorized modification, lack of maintenance or deliberate damage to the machine is not allowed by law. And this could inhibit the warranty.

Autoclaves cannot be used for fluids that can produce vapours which can form inflammable or explosive mixtures.

Use it only for the intended use of the equipment as specified in the catalogue, or manuals.

2 Warranty period

This autoclave is guaranteed for one year.

Warranty does not cover damages caused by an incorrect use or causes beyond the control of J. P. SELECTA, s.a.u.

Any manipulation of the equipment by personnel not authorized by J.P. SELECTA, s.a.u. automatically cancels the guarantee.

3 Technical Specifications

Autoclave volume:	50L	80L	150L
Mains Volt:	230V 50/Hz		400V III/50Hz
Electrical power: (W)	3000	4850	7500
Max. Current: (A)	13	21	20
Fuse: (Fast. Type H, dimensions)	20A	25A	32A
Mains wire	PHASE Neutral, Ground (sec.2,5mm ²) R,S,T neutral, earth (sec 6mm ²)		
Weight (Kg Net):	108	136	250
Vessel Dimension usable (cm Ø x high)	Ø30x63	Ø40x60	Ø50x70
External Dimensions: (cm DepthxWidthxHeight)	117x48x62	119x58x72	118x80x90
Vessel capacity: (L)	50L	80L	150L
Used volume:	48L	78L	140L
Max load: (Kg) (Metal)	45		50
Weight per support (Kg/cm ²) Empty	25		
Weight per support (Kg/cm ²) Load + Water	35		
Total space + doors open	Add 1m at height		
Tap Water in max. pressure allowable (bar)	2		
Tap water in min. pressure allowable (bar)	0,5		
Water vol. per cycle (L)	0,75	1,50	5.0
Water quality suggested	Conductivity demineralized > 100 microS/cm		
Thermal heat to ambient: (W/h)	750	850	850
Water reservoir volume: (Litres)	10		15
Water in connector	3/4"		
Air entry filter porosity	0.20 micres		
Noise level: (1m)	Lower than 70 dBA		
Operating time:	Continuous left 20 minutes between cycles.		
Ambient:	Indoor use		
	High up to 2000m		
	Ambient operating range from 5°C to 40°C.		
	Maximum relative humidity 80% for temperatures up to 31°C, linearly decreasing to 50% relative humidity at 40°C		
Over voltage level:	Category II		
Contamination degree:	2		
Steam contact material	Stainless Steel, Copper, Teflon, Brass		
Vessel material:	Stainless Steel		
Frame material:	Coated steel		
Steam generator:	Stainless Steel		

4 Packing List

Autoclave delivered together with:

- Instruction manual
- Internal rack of bottom of the tank

5 Overview and recommendations

These autoclaves can sterilize solids, with or without being wrapped, liquids, culture media and also have disinfecting drying cycles.

The material must be clean, free of any kind of scale or residue.

Do not overload the trays, racks or baskets and leave a space between them of 1 or 2 cm, in order to allow the passage of steam and facilitate drying.

5.1 Purge

The purging operation is performed by the autoclave to remove air inside the boiler as so fill the vapour space.

The purge system must be suitable to sterilize intended load.

5.2 Atmospheric purge

In atmospheric purging, air comes out gravimetrically, impulsed by steam during a certain time.

Suitable for solid unwrapped (Tweezers, needles, scalpels ...).

5.3 Pulsed vacuum purge

A vacuum pump is used to extract air from the boiler. During purging, one or more extractions purged air is produced. This system is much more effective than atmospheric purging.

Suitable for wrapped loads, catheters, ...

5.4 Bowie Dick test and Vacuum Test

These tests are pre-programmed cycles.

Bowie Dick test used to verify operation of the purging by vacuum fractional.

The Vacuum Test is designed to verify the operation of the vacuum pump and sealing the autoclave. This test must be carried out with a cold equipment (it is usual to do it before starting the work sequence).

5.5 Sterilization of unwrapped material

You have to place the instruments on waterproof paper and they do not have to touch each other.

5.6 Sterilization of wrapped materials and pipes

You only have to place one instrument per bag.

Placing the bags upright they do not touch one another nor the walls of the autoclave.

The tubes must be clean, rinsed and drained. Ensure that the tubes do not form folds that obstruct the inner diameter.

5.7 Sterilization of vessels

Never place sealed containers.

Place the containers upside down to prevent water tanks.

5.8 Sterilization of liquids

Place fluids in suitable containers to work at sterilization temperature and place them above trays so that they can collect spills.

The containers were filled about 2/3 full. Do not exceed this value. Do not close them tightly, they will be plugged with cotton or some sort of cap to facilitate exhaust of the pressure vessel and not to form pressure inside them.

Avoid using containers with narrow necks.

At the end of the sterilization cycle, the autoclave was let cool freely (until pressure is 0 kg / cm²).

5.9 Water to use

Type of water is of vital importance in the operation of the equipment and its subsequent maintenance.

It is recommended to use distilled water with a conductivity level higher than 100 uS / cm so that the steam generator water level sensor works correctly.

Mains water may contain calcareous and ferritic residues that can damage the passage through valves and heating components. Use soft water (see table) to avoid encrustations.

It is also not advisable to use pure or ultra pure water because it is very corrosive and can even damage the stainless steel of the chamber.

Water type	mg/l	*fH	*dH	*eH
Soft water	≤17	≤1.7	≤0.95	≤1.19
Slightly hard water	≤60	≤6.0	≤3.35	≤4.20
Moderately hard water	≤120	≤12.0	≤6.70	≤8.39
Hard water	≤180	≤18.0	≤10.05	≤12.59
Very hard water	>180	>18.0	>10.05	>12.59
mg/l: milligrams per liter of calcium carbonate water (CaCO ₃) *fH: French grade (10.0 mg CaCO ₃ /l) *dH: German grade (17.8 mg CaCO ₃ /l) *eH: English grade (14.3 mg CaCO ₃ /l)				

5.10 Preheating

When starting the equipment, to accelerate the sterilization cycle, it heats the internal walls of the chamber. The temperature parameter can be selected from the SETTINGS menu.


6 Equipment description

6.1 Control panel elements (50 & 80L)


1. Touch screen
2. Opening and closing lid lever
3. Optional printer accommodation
4. Main power switch

18. Clean water drain valve


5. Door joint
6. Water tank filling hole
7. Access port inside the boiler


19. Temperature probe
20. Mesh filter drain


6.2 Rear part elements (50 & 80L)

8. External water source inlet
9. Safety valve
10. Steam manual valve
11. Mains power socket
12. Aeration filter
13. Drainage selection valve
14. Condensed steam outlet
15. Liquids mode outlet
16. RS232 connector
17. USB connector



6.3 Front part elements (150L)

1. USB connector for data recorder.
2. Touch screen.
3. Lid lock.

4. Mains switch.

16. Mesh filter drain.

5. Water tank refilling entry.
6. Gasket.


6.4 Rear part elements (150L)

7. Connector RS232 (Serial port)

8. Safety valve
9. Steam output manual valve
10. Water tap entry
11. Aeration filter
12. Maximum thermostats
13. Condensate outlet
14. Tank water outlet
15. Drain valve

6.5 Cycle types and temperature range

The autoclave allows you to configure the sterilization cycle according to the characteristics of the load to be sterilized:

Program type	Purge type	With internal probe	With external probe	Material type
SOL1	Atmospheric	105°C-134°C	105°C-130°C	Unwrapped instruments
SOL2	Vacuum 2 fractions	115°C-134°C		Bagged instruments
SOL3	Vacuum 3 fractions	115°C-134°C		Double bagging, porous
LIQ	Atmospheric	105°C-134°C	105°C-130°C	Liquids, culture media

Sterilizing time range:

From 3 to 60 minutes

Drying time range:

From 0 to 60 minutes

Temperature resolution:

±1°C

Time resolution:

1 minute

6.6 Pre-set cycles

Autoclave incorporates 9 pre-set working cycles which can be modified according to the user needs.

Program No.		Temperature °C	Time (minutes)	Dry	Purge
0	Not wrapped load	134	4	10	2x vacuum cycles
1	Wrapped load	134	12	30	3x vacuum cycles
2	Prions	134	20	20	3x vacuum cycles
3	Delicate not wrapped	121	15	15	2x vacuum cycles
4	Delicate wrapped	121	30	30	3x vacuum cycles
5	Instrumental	126	12	25	3x vacuum cycles
6	Disinfection	105	25	15	Atmospheric
7	Liquids	121	30	-	Atmospheric
8	Solids	121	15	30	Atmospheric
9	Cleaning program*	134	4	30	Atmospheric

* Fixed program, it can not be modified.

6.7 Test cycles

Also, the autoclave Autester ST DRY PV III incorporates two test programs.

Bowie Dick test	134°C	3' 30"
Vacuum Test	0,7 bar	5+ 10'

7 Installation

Place the autoclave near a power outlet suitable for the machine.

The autoclave will sit on a flat, horizontal, stable surface, leaving a space of 30 cm around the machine.

Immobilize the autoclave, adjusting it to the ground, by using the stops located on the front of the equipment, provided for this purpose.

Open the lid. Place the rack at the bottom of the boiler. Check that there are no packaging elements inside the boiler. Verify that the boiler is clean and it does not contain particles that could damage the passage through the valves.

Attach a tube in the "OUT CONDENSED" spike to collect condensate formed during the phases in which the vacuum pump works. Do not place the vessel at a higher level of the nipple.

If you plan to work in Liquid mode, attach a tube to the "OUT LIQUIDS" spike. Fill the auxiliary tank with water.

7.1 External water inlet.

The autoclave can be filled through the external network through a 3/4 "socket as seen in the picture. **The supply pressure must never be higher than 2bar.**

This form of filling must be set from the configuration panel.

7.2 Selection of drainage system

Place the valve control pointing it towards inside the machine with the selector valve indicated, **if you want to recover the water at the end of each sterilization process.**

However, place the valve control in the direction of outside, **if you want to expel water outside**, at the end of each sterilization process. Fix a hose on the valve outlet nipple and also fix the other end to a container or drain, taking care **NOT TO OBSTRUCT THE TUBE.**

ATTENTION

Steam outlet in case of overpressure.

Risk of burns.

Position this equipment in the proper way to prevent this risk from occurring.

Autester ST DRY III 50-80L

Autester ST DRY III 150L

Set **firmly** drain hoses on these three outputs !!




ATTENTION

Do not use the external drain option without connecting first a pressure hose. Fix the hose so that it cannot move due to the pressure at the end of the cycle. Do not hold the hose with your hands, you may suffer burns.

7.3 Measurements of the external autoclave connections

8 Operation

Plug the cable into the base of the power outlet.


The autoclave lid only opens with the autoclave plugged into the network (illuminated switch)

Do not FORCE the opening lever

Switch on the device by activating the main switch (4) and the main screen will appear.

8.1 Main menu

All the functions of the control panel are selected by the touch screen. Press with your finger or with a pointer in the center of the touch screen.

AJUSTES
SETTINGS

PROGRAMACIÓN
PROGRAMMING

TRABAJO
TASK

LANGUAGE DATE ENTRY PASSWORD

PRINTER UNITS CHANGE PASSWORD

PROBE WATER

SETUP HH:MM:SS DD/MM/YY

SELECT LANGUAGE:

ENGLISH  

SETUP > LANGUAGE HH:MM:SS DD/MM/YY

POINT TO SET TIME / DATE

TIME 15 : 01

DATE 11 / 05 / 15

SETUP > DATE / TIME HH:MM:SS DD/MM/YY

ACTIVATE PRINTER NO  

PRINTER INTERVAL 1 min  

PRINTER FEED

SETUP > DATE / TIME HH:MM:SS DD/MM/YY

PROBE

INTERNAL  

SETUP > PROBE HH:MM:SS DD/MM/YY

TEMPERATURE UNITS.

°C  

SETUP > UNITS HH:MM:SS DD/MM/YY

8.2 Parameters settings

To set operation parameters as: language, temperature units, control probe, etc...

Point over «SETTINGS».

Point over the parameter to be modified.

8.2.1 Language:

Set language between:

Spanish.
English.
French.

8.2.2 Date and time:

Press over icons to modify day, month, etc...

Press  to return to the main menu or press  to modify other parameters.

8.2.3 Printer:

(Optional accessory, must be factory assembled)
Enables / disables the printer and the data writing time.
"Printer Feed" feeds the printer paper.

Press  to return to the main menu or press  to modify other parameters.

8.2.4 Internal or external (heart) probe selection:

(Optional accessory, must be factory assembled)

Internal probe: The autoclave controller takes the vessel water temperature of the autoclave's bottom to control the cycle.






External probe: (accessory). Also called «heart probe», it takes as reference the temperature of the interior of the load during the sterilization phase.

Press  to return to the main menu or press  to modify other parameters.

8.2.5 Temperature units:

Choose units between °C or °F.

Press  to return to the main menu or press  to modify other parameters.

8.2.6 Activates / disables external water entry

With this parameter set in "**activate**" the autoclave, the level of the water tank for steam generation will always be maintained in working limits. Manual filling will not be necessary.

With the "**disable**" value set, the autoclave takes water for generating steam from the autoclave tank itself. Filling additional deposit will be manual.

Press to return to the main screen or to modify another parameter.

8.2.7 Parameters protected by a password

The password protects the parameters that may affect the operation of the autoclave.

To introduce the correct password, press over «INSERT PASSWORD».

- **Calibration of measuring devices (password "0000"):**

Internal temp: offset value (Offset) of the water probe

External Temp: offset value (Offset) of the heart probe

Pressure: displacement value (Offset) of the pressure probe

Pre-heat: preheating temperature of the boiler

Program 0-8 Temperature (default 95°C)

Log interval: time between the registration of two parameters 1-60min

8.3 Test cycles (Vacuum test & Bowie Dick)

The autoclave has two EN13060 test cycles programs to made easy operation qualification:

- Vacuum test: Qualifies the autoclave tightness and its capability to get suitable vacuum according to EN13060 (This test must be carried out with the cold boiler)
- BOWIE DICK test: Qualifies the steam quality and steam capability to penetrate into wrapped loads.

This test cycles should be run periodically, and results must be stored as quality system records.

8.4 Sterilization programs

Up to 9 cycles could be stored on the autoclave, editable by the user.

The cycle sets the treatment to apply over samples, by means of the setting parameters.

The 9 programs parameters are loaded by default in the device. See the table in section 6.4 for the default values.


8.4.1 Modify & save a program

- Select the cycle to modify by pressing over:

- The cycle (program) sets the sterilization treatment applied to the load, by means of four parameters.

- Purge type.
- Sterilization temperature.
- Sterilization time.
- Drying time.

Note: Checking this box which is the operating program when starting the autoclave. Select another program and press the SEL key to select another of the programs in the list.


Purge type

• Setting purge type

Press on the purge type to select:


Solid 1: Atmospheric purging	For unwrapped loads.
Solid 2: 2 vacuum cycles.	For wrapped loads.
Solid 3: 3 vacuum cycles.	For double wrapped loads, catheters
Liquids: Atmospheric purging	Culture media in demi-close bottles, without drying.

Sterilization temperature

• Setting sterilization temperature


Press over temperature label to set between 105°C ... 134°C.
(In the case of working with a heart probe, the maximum temperature limit is 130°C)

Sterilization time

• Setting sterilization time

Press over sterilization time to set time from 3 to 60 minutes.


Drying time

• Setting drying time


Press over drying time to set time from 0 to 60 minutes.
(0 min = without drying)

(If the purge type for the selected program is LIQ, this option disappears, it is incompatible with the purge type)

Save program changes

• Store programs

Quit to main menu

• Return to main menu

Note: Checking this box  which is the operating program when starting the autoclave. Select another program and press the SEL key to select another of the programs in the list.

SETTINGS	PROGRAMMING
TASK	

8.5 Load

Open the lid. Check tank water level.

Load the previously cleaned and settled material baskets.

Check that the seal (boiler top) is in place.


If it is installed and if it is necessary, place the heart probe in the chosen reference point.

Close the lid.

8.6 Select Program


Select the suitable program for sterilization and drying the load.

See chapter 8.4 to set program cycle.

Drying time

Shows running
program status


8.6.1 Application of the cycle

During the cycle (also called program) autoclave screen will indicate the phase that is being made.


Phases are:

WATER:	Water inlet to the boiler.
VACUUM:	Air extraction.
HEAT:	Steam generation.
STERILE:	Sterilization of the load.
EMPTY:	Draining the boiler steam and water.
UNPRESS:	Depressurization boiler.
OK:	Sterilization conditions have been achieved according to cycle parameters.
FAIL:	Not achieved all the conditions for the sterilization.


This is the "End of cycle" message.

8.6.2 Run a pre-programmed cycle


1 Press over «PROGRAMMING».


2 Select a cycle with


PROGRAM ... appears and also the memorized parameters in the displayed program.


3 Once selected, press on and confirm.

4 Press over to go to the main menu.

5 Press over «TASK». The program we have chosen appears. Check the door is closed.

6 Press over «RUN» to start.

```
*****
14/05/2015 09:52
=====
PROG: 0
MODE: SOLID2
TEMP: 134°C
TIME: 04'
DRY : 10'
=====
LOT: 000010
*****
```

TIME TEMP. PRESS. STATUS

```
023°C +0.0bar VACUUM
021°C -0.2bar VACUUM
021°C -0.4bar VACUUM
021°C -0.5bar VACUUM
022°C -0.5bar HEAT
023°C -0.5bar HEAT
023°C -0.5bar HEAT
025°C -0.5bar HEAT
060°C -0.3bar HEAT
071°C -0.1bar HEAT
080°C +0.0bar HEAT
089°C +0.2bar HEAT
095°C +0.3bar HEAT
100°C +0.5bar HEAT
095°C +0.0bar VACUUM
090°C -0.1bar VACUUM
086°C -0.3bar VACUUM
083°C -0.4bar VACUUM
081°C -0.4bar VACUUM
078°C -0.5bar VACUUM
076°C -0.5bar VACUUM
079°C -0.4bar HEAT
088°C -0.2bar HEAT
095°C +0.0bar HEAT
101°C +0.1bar HEAT
106°C +0.3bar HEAT
110°C +0.5bar HEAT
114°C +0.7bar HEAT
119°C +0.9bar HEAT
122°C +1.2bar HEAT
126°C +1.4bar HEAT
129°C +1.7bar HEAT
133°C +2.0bar HEAT
134°C +2.0bar HEAT
03' 134°C +2.0bar STERIL
02' 134°C +2.1bar STERIL
01' 134°C +2.1bar STERIL
00' 134°C +2.1bar STERIL
132°C +1.9bar EMPTY
124°C +1.1bar EMPTY
112°C +0.4bar EMPTY
104°C +0.1bar EMPTY
00' 099°C +0.0bar DRY
01' 093°C -0.2bar DRY
02' 089°C -0.3bar DRY
03' 085°C -0.4bar DRY
04' 082°C -0.4bar DRY
05' 079°C -0.5bar DRY
06' 076°C -0.5bar DRY
07' 073°C -0.6bar DRY
08' 071°C -0.6bar DRY
09' 069°C -0.6bar DRY
068°C +0.0bar UNPRES
RESULT: PASS
```

9 USB data recorder

The equipment has an internal memory for data acquisition (Datalogger). Once the cycle is finished, it is possible to download this information to a USB memory through a USB connector.

The device (USB pen drive) must be formatted to FAT 32.

Operation

- By default, the logger is active with a 1 minute logging interval. Follow the instructions to modify the recording interval (see 8.2.7).

With the cycle completed:

- Connect the «USB pen drive» to connector.

- When you connect a "USB pen drive", the USB symbol will appear in the lower bar of the screen.

- From the touch screen select:
SETTINGS> DOWNLOAD LOGS

A download button will appear and by pressing it all logs will be downloaded and the internal memory emptied.

- The screen will indicate when the recording of the logs has finished. Disconnect the «USB pen drive» from the device.

- Connect the «USB pen drive» to a computer to download the data file.

Important Note: The «USB pen drive» must have a device name. Can not be blank.

The created file is of text type, that can be opened with WordPad, WORD or EXCEL.

Note: If the data is not downloaded, when the device memory is at 80% capacity, a warning message will appear indicating that the memory needs to be downloaded.

If the message is ignored and the memory is completely filled, the oldest cycles will be erased.

Approximately the memory has capacity to store about 150 cycles (it will depend on the total time that each cycle has).

Example of registration obtained by the printer

ERROR 10
WATER RESERVOIR EMPTY
wait until pressure
and autoclave temperature
becomes low and press OK

OK

ACTUAL: +0.00BAR 23°C HH:MM:SS DD/MM/YY

10 Alarms

If a failure occurs during operation, a message will appear. These failures are indicated by showing ERROR on the screen.

Please contact Technical Service or a distributor if you get one of them.

ERROR	MESSAGE	DESCRIPTION
01	E2PROM FAIL	EEPROM faulty error. Stop and restart the equipment. If the error persists, contact an official technical service.
02	PRESSURE SENSOR FAIL	Pressure sensor error. Failure sensor or disconnected cable.
10	WATER TANK EMPTY	Water reservoir empty. Fill the equipment water tank.
12	DOOR OPEN	Door open. Close the door to start the computer.
13	OVERPRESSURE	The pressure has exceeded maximum admissible pressure of +2.5bar .
14	VACUUM FAIL	Vacuum error. Check the manual depressurization valve is in working position. Vacuum leakage or pump failure.
15	PRESSURE FAIL	Pressure error during purge phase. Insufficient water or water heater failure.
16	INTERNAL OVERTEMP	Over temperature error. Internal temperature has exceeded maximum limit of 150°C .
17	EXTERNAL OVERTEMP	External heating element overtemperature. External temperature has exceeded maximum limit of 150°C .
18	STERILIZATION FAIL	Sterilization error. Set point margins overcome. There is air inside.
19	SETPOINT UNREACHABLE	Set value error not reached. No water in the system or heater failure.
20	EMPTY TOO LONG	Empty error, takes too long to empty the system. Valves or drain blocked or failure.
21	PRINTER PAPER END	No paper error, printer requires paper.
22	PRINTER ERROR	Printer error. Possible paper jam.
23	POWER FAIL	Power error. Detected in the start-up. The equipment has been switch off before the end of the cycle.
26	TEST FAILED	Vacuum test error.
29	SAFETY THERMOSTAT	Thermostat open. Overheating. Unlock thermostat and check water level
30	LEVEL ERROR	Steam generator sensor error or blocked tubes. More than 200s filling water.
31	PURGUE ERROR	Maximum purged time exceeded error.
32	DEPRESS FAIL	Maximum depressurization time exceeded error.
33	EXT. PROBE FAIL	Second temperature probe failure in LIQUID mode.
34	INTERN UNDERTEMP	Low temperature. Internal temperature does not reach 8°C .
35	COMMUNICATION FAIL	Communication error between touch screen and autoclave. This error message will appear when, after several attempts, messages do not receive any reply.

11 Safety Systems

11.1 Safety pressure switch

When an abnormal over pressure is produced in the sterilization chamber, the safety pressure switch acts in a way to stop the autoclave and it opens the aeration electrovalve.

11.2 Safety valve

If an abnormal over pressure is produced, the safety valve is opened, expelling the steam.

CAUTION: the steam valve must be oriented in a direction that can not reach the personnel.

This valve is calibrated to initiate the opening between 2.2 kg/cm² and 2.5 kg/cm².

In case of safety valve reset, the machine must be prepared to advise technical personnel.

11.3 Safety thermostat


When there is an abnormal temperature rise or a significant decrease in the water level (with consequent overtemperature in the steam generator), the safety thermostat will start to stop heating.

In the 150L model, the thermostat must be reset to restart the cycle. To reset the thermostat, press the white button, then a click will be heard (see attached photo). A message will appear on the screen.

11.4 Unblocking the door

In case of open lid blocked by a safety system failure make as follows:

- 1 Stop the machine.
- 2 Disconnect the machine from the mains.
- 3 Place an outlet and a tube in the aeration valve.
- 4 Carefully open the valve to avoid steam leakages. Be careful to avoid possible burns.
- 5 Wait for the total steam emptying.
- 6 Extract the upper cover where display is located in order to put it aside (it is not necessary to disconnect any wire). It does not use screws, pull tightly.
- 7 Remove the plastic cap and pull slightly until the lever is released.

12 Maintenance

Before removing the lid, disconnect the machine from the mains.

The control board should only be changed by a technical service or authorized personnel.

CLEANING:

- Wall panels and stainless steel: Use steel cleaning home products.
- Limescale deposits: Use home cleaning products.
- TFT Screen: use a cloth and a glass-cleaning liquid
- Plastic lids: use a cloth and a glass-cleaning liquid
- Perform a cleaning program (No. 9) - See ***Note**

Monthly operations:

- Clean front bottom drainage filter.
- Change re-used water from water reservoir. Empty the water tank with valve 18 (in the 50 & 80L autoclaves)/16 (in the 150L autoclave)
- Record results from vacuum test and Bowie Dick test.

Three months operations:

- Change air entry filter (autoclave rear side)

Yearly operations:

- Check and contrast temperature autoclave against traceable temperature standard.
- Check safety valve operation.
- ATTENTION: this operation should be run for qualified personnel. Dangerous contact with hot steam.
- Visually check the state of electrical wiring, ground and mains.






***Note** - Cleaning program:


- Make sure the tank is empty.
- Select the LIM program as indicated in point 8.6.2.
- Fill with water manually (about 2 liters of water) and pour through the tank door.
- Close the lid.
- Run the program normally.

It is recommended to run this program after a heavy work day or once a week to avoid deposits in the pipes.

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
**Notice to customers about recycling:**

The product is made up of various components and various materials that must be recycled or, failing that, deposited in the corresponding debris removal sites when the product's life has been completed or when otherwise it is necessary to dispose of it. To do this, the end user who acquires the product must know the current regulations of each municipality and / or locality based on the waste electrical and electronic equipment. The user who acquires this product must be aware of and responsible for the potential effects of the components on the environment and human health as a result of the presence of hazardous substances. Never place the product in a conventional container of citizen scope if a previous dismantling and knowledge of the components that incorporates. If you do not know the procedure to follow, consult with the city council for more information.