




▶

KBN - Weighting accesories

▶ Watch glass ground edge

- 1 |
- These Endo Glassware watch glasses are concave, made of glass pieces that present ground edges and are available in different sizes.
- 2 |
- They can be used as weighing trays, for drying wet solids, for the observation of small invertebrates under a stereomicroscope, as recipients in oil tasting or as a lid for beakers, flasks and other recipients.



Code	Diameter (mm)
KBN001	45
KBN002	50
KBN003	60
KBN004	70
KBN005	80
KBN006	90
KBN007	100
KBN008	120
KBN009	150
KBN010	180

▶ Watch glass burnt edge

- 1 |
- Watch glasses with burnt edges to avoid cuts. Specially recommended to evaporate liquids or weigh samples, they are available in different sizes.



Code	Diameter (mm)
KBN011	40
KBN012	50
KBN013	60
KBN014	70
KBN015	80
KBN016	100
KBN017	120
KBN018	140
KBN019	160