



Analog multi-position, RSLAB-9

Code LJG002, LJG001

- Allows for a gentle yet effective mixing of biological samples. Two different models are available: the first model is disk-shaped, and the second model resembles a rotisserie.
- Compact design enables the equipment to be placed inside an incubator or refrigerator.
- Compatible with a wide variety of accessories.
- Protection according to DIN EN60529: IP21.
- Supplied with:
 - 3 Bags of 12 adapters for 1.5 mL tubes.
 - 1 Bag of 12 adapters for 15 mL tubes.
 - 1 Bag of 12 adapters for 50 mL tubes.



Code LJG002



Code LJG001

Technical specifications

| Code | LJG002 | LJG001 |
|-------------------|----------------|--------------------|
| Model | RSLAB-9 Disc | RSLAB-9 Rotisserie |
| Speed range | 0 – 80 rpm | 0 – 80 rpm |
| Shaking movement | Rotational | Rotational |
| Operation type | Continuous | Continuous |
| Aprox. dimensions | 300x220x310 mm | 530x150x190 mm |
| Aprox. weight | 2,0 Kg | 1,7 Kg |
| * Brand: DLAB | | |

Accessories included with the equipment



TUBE CLAMPS

| Code | Description | Pack |
|--------|---------------------|------|
| LJJ011 | For tubes of 1,5 mL | 12 |
| LJJ012 | For tubes of 15 mL | 12 |
| LJJ013 | For tubes of 50 mL | 12 |



Digital multi-position, RSLAB-9

Code LJG004, LJG003

- Digital control of speed and time on the LCD screen. Smooth and effective rotational movement, equipped with a braking device.
- Wide range of accessories for using different types of containers, ensuring stable fixation during agitation.
- Protection according to DIN EN60529: IP21.



Code LJG003



Code LJG004

Technical specifications

| Code | LJG004 | LJG003 |
|-------------------|---------------------------|----------------------------|
| Model | RSLAB-9 digital disc | RSLAB-9 digital rotisserie |
| Display | LCD | |
| Speed range | 10-70 rpm | |
| Shaking movement | Rotational | |
| Shaking angle | Adjustable from 0° to 90° | |
| Operation type | Continuous or scheduled | |
| Timer | 1 - 1199 min | |
| Aprox. dimensions | 510 x 190 x 260 mm | |
| Aprox. weight. | 8,0 Kg | |




Digital multi-position, RSLAB-9

Code LJJ004, LJJ003

Accessories included with the equipment



Rotisserie-type, vertical
Code LJJ006*
For 16 tubes of 50 mL
**included with the equipment*



Disk-type
Code LJJ009*
For 8 tubes of 50 mL
** included with the equipment*

Accessories compatible with the equipment



Rotisserie-type, horizontal
Code LJJ001
For 48 tubes of 1,5 mL



Rotisserie-type, horizontal
Code LJJ002
For 24 tubes of 15 mL



Rotisserie-type, horizontal
Code LJJ003
For 24 tubes of 50 mL




Rotisserie-type, vertical
Code LJJ004
For 32 tubes of 1,5 mL



Rotisserie-type, vertical
Code LJJ005
For 16 tubes of 15 mL



Disk-type
Code LJJ007
For 60 tubes of 1,5 mL



Disk-type
Code LJJ008
For 16 tubes of 15 mL



4 support rods
Code LJJ010
To integrate another disk-type
support