



ENMIND EN-V5 Vet Infusion Pump

Code ZMZ015

AXAVET's infusion pump allows precise and safe administration of liquids and medicines.
Quick access to the 20 most recent configurations.

- Small and lightweight
- Large touch screen for clear display of parameters including remaining infusion time
- Different settings available, such as: total perfusion volume, required flow rate, total perfusion time
- High sensitivity of settings (bubbles and occlusion pressure)
- High-performance battery (up to 10 hours of operation)
- Safety and operational alarms
- Self-diagnosis function
- Suitable for transfusions
- Compatible with standard droppers. Holder for serum carrier included



Technical parameters

Code	ZMZ015
Accumulated volume	From 0,1 to 9999 mL
Touch screen	4,3 -inch color LCD
Flow rate range	0,1 – 2000 mL/h
Modes:	
Flow rate and volume	
Total time and volume	
Time and total flow rate	
Body weight	
Stored modes	
Micro mode	
Accuracy	±5%
Air bubble detection	20, 50, 100, 200, 400, 800 µl
Occlusion settings	5 levels (50, 150, 300, 600, 900 mm Hg)
Data review	5000 entries
Protection class	IP44
Power supply	100-240VAC, 50/60 Hz
Dimensions (LxWxH)	14x13x9 cm
Weight	1,45 kg