

COMECTA Stainless steel receptacles.

QUALITY AISI 304.

LABORATORY-HOSPITAL + SERIES



ROUND BASINS

Part No.	Ø exterior	Ø interior	Height	Volume
8081810	133 mm.	118 mm.	45 mm.	400 ml.
8081812	175 mm.	157 mm.	58 mm.	850 ml.
8081814	217 mm.	199 mm.	66 mm.	1650 ml.
8081815	255 mm.	245 mm.	110 mm.	3400 ml.
8081817	318 mm.	297 mm.	117 mm.	7000 ml.



STERILIZATION BOXES WITH LID

Part No.	Height	Width	Depth
8085504	30 mm.	170 mm.	70 mm.
8085505	50 mm.	200 mm.	100 mm.



RECTANGULAR TRAYS

Conical shape

Part No.	Exterior Height	Exterior Width	Exterior Depth	Interior Width	Interior Depth
8080801	40 mm.	190 mm.	125 mm.	158 mm.	112 mm.
8080802	45 mm.	244 mm.	161 mm.	211 mm.	147 mm.
8080803	50 mm.	305 mm.	195 mm.	270 mm.	177 mm.
8080804	55 mm.	338 mm.	261 mm.	301 mm.	241 mm.

Long shape

Part No.	Exterior Height	Exterior Width	Exterior Depth	Interior Width	Interior Depth
8080208	18 mm.	175 mm.	62 mm.	145 mm.	49 mm.
8080207	28 mm.	380 mm.	160 mm.	351 mm.	145 mm.



BOXES WITH ROUND LIDS

Part No.	
8085514	140 Ø x 65 mm high. Volume 1000 ml.
8085515	180 Ø x 75 mm high. Volume 1500 ml.



KIDNEY DISHES

Part No.	
8085301	32 high x 200 mm long.
8085302	37 high x 250 mm long.
8085303	42 high x 300 mm long.



TALL BEAKERS AND JUGS

Mirror finish stamped beakers and jugs with a curved lip.
Easily cleaned and resistant to large numbers of chemical agents.

1 Standard	2 With Spout	3 With handle and spout		
Part No.	Part No.	Part No.	Ø x Height	Volume
8002902	8012902	8092902	72 x 84 mm.	250 ml.
8002905	8012905	8092905	86 x 129 mm.	500 ml.
8029010	8129010	8929010	102 x 150 mm.	1000 ml.
8029020	8129020	8929020	132 x 184 mm.	2000 ml.
8029030	8129030	8929030	160 x 180 mm.	3000 ml.
8029040	8129040	8929040	180 x 200 mm.	4000 ml.
8029050	8129050	8929050	200 x 220 mm.	5000 ml.



LIDS FOR BEAKERS AND JUGS

Part No.	Beaker size	Part No.	Beaker size
8000250	250 ml.	8003000	3000 ml.
8000500	500 ml.	8004000	4000 ml.
8001000	1000 ml.	8005000	5000 ml.
8002000	2000 ml.		



BUCKET WITH HANDLE

Buckets, Graduated from 1 to 12 litres.
320 Ø x 300 mm high. Volume 12 litres.

Part No. **8000252**

Bucket lids

Part No. **1002521**




MEDICINE CUPS

Part No.

8085530 Graduated, Capacity 30 ml.

8085531 Graduated, Capacity 50 ml.

8085532 With handle, no graduations. Capacity 300 ml.



SERILISATION JAR

100 Ø x 155 mm high.
Volume 1000 ml.

Part No. **8080904**



ROUND BASIN

370 Ø x 145 mm high. Volume 9000 ml.

Part No. **8000140**



SUPPORT FOR ROUND BASIN

Portable with rubber wheels and brake.

Basin Part No. 8000140.

Part No. **6010140**



SUPPORT FOR BUCKETS

Portable with rubber wheels and brake.

For bucket Part No. 8000252.

Part No. **6010252**