



items for clinical use  
[Z] OTHERS



## ZMB035 - Hemacytometer



1| Nahita hemacytometer is made up of a complete set with all the necessary elements to determine erythrocyte and leukocyte density in blood samples.

2| Its counting chamber is one glass piece construction, ensuring long term durability and resistance.

3| Furthermore, it presents a double 0,1 mm-deep Neubauer ruling for cell counting and it is also equipped with two special 0,4 mm-thick cover glasses.

4| The set also includes two blood diluting pipettes (for erythrocytes and leukocytes), each one with its corresponding rubber tube of approximately 35 cm of length, and its colour-coded mouthpiece for an easier identification.

5| The complete set is supplied in a plastic case to keep it away from dust and perfectly stored.

6| Without IVD certification, for sales outside the EC

ITEMS FOR CLINICAL USE  
[Z] OTHERS

Note: manufacturer reserves the right to modify features and photographs of the goods