



WIRE BASKETS


Made of AISI 304 stainless steel sheet with electro-bright finish. Designed to fit the Selecta autoclaves. (See page 53).

Part No.

1004771 220 Ø x 140 mm high.
1004772 220 Ø x 240 mm high.
1000495 265 Ø x 185 mm high.

Part No.

1000496 365 Ø x 265 mm high.
1000780 435 Ø x 205 mm high.



ACCESSORY


Basket dividers, stainless steel AISI304, electrobright finish.

Part No.

1001217 4 compartment divider
1001218 6 compartment divider
1001222 7 compartment divider

For baskets, Part No.

1000495
1000496
1000780



WIRE BASKETS

Electropolished stainless steel AISI 304.

Part No.

Dimensions

1025115 210 Ø x 180 mm high
1025121 150 Ø x 120 mm high
1025124 240 Ø x 180 mm high

With cross dividers.

Part No.

1125115
1125121
1125124



STERILIZATION DRUMS

Made from AISI 304 stainless steel. Designed to fit the Selecta autoclaves. (See page 53.)

Part No.

Ø cm

Height cm

1002418 28 18
1002419 35 27
1002420 28 18

Doesn't match Selecta autoclaves.



STERILIZER CAN FOR PIPETTES

Made from Stainless steel AISI 304.

Part No.

height

Ø

1000829 400 72
1001126 400 130

STERILIZATION CONTAINER FOR PETRI DISHES

Made of AISI 304 stainless steel. Complete with Petri dish rack for up to 120 mm Ø.

Models

height / Ø

No. of

Part No. mm plates
1001127 260 / 130 10
1001128 400 / 130 18



PETRI DISH RACK

Made of AISI 304 stainless steel, with electro-bright finish.

For 10 dishes of 80 to 100 mm Ø.

Part No. **1000259**



PIPETTE BASKET

Made of AISI 304 stainless steel with electro-bright finish. 90 mm Ø x 250 mm high.

Part No. **1000251**



STERILIZATION CONTAINER FOR PIPETTES

Made from Stainless steel AISI 304. Lid interior and base lined with silicon to protect pipette end and tip.

Part No.

height mm

Ø mm

1001597 290 65
1001598 340 65