


**OPTIONS:** for vertical autoclaves, need to be factory fitted.




**Pt 100 Sensor.** 2<sup>nd</sup> Probe for reading sample temperature.

To guarantee in all circumstances that the system has completed the cycle to the set time and temperature parameters, the Autester ST DRY PV III autoclave can be installed with a temperature probe that can be placed within the culture flask, residual packets etc.

When the probe is connected to the controller of the autoclave, it must be selected in settings, whereby the real time temperature values detected by the probe, are displayed and printed.

For autoclaves PRESOCLAVE III and PRESOCLAVE PLUS this Pt100 sensor only reads the temperature with an external thermometer. Part No. **4001219**

## ACCESSORIES FOR VERTICAL AUTOCLAVES



AISI 304 stainless steel wire baskets.

Part No.	Ø	Height cm	Suitable for autoclaves Part No.	Number of baskets
<b>1004771</b>	22	14	<b>4001745</b>	1
<b>1004772</b>	22	24	<b>4001756</b>	1
<b>1004773</b>	23	17	<b>4001757</b>	2
<b>1000495</b>	27	20	<b>4001758 - 4001770 - 4002516</b>	3
<b>1000496</b>	36	28	<b>4001759 - 4001771 - 4002519</b>	2
<b>1000780</b>	44	22	<b>4001772 - 4002518</b>	3
<b>1000490*</b>	26	19	<b>4001758 - 4001770 - 4002516</b>	3
<b>1000491*</b>	36	26	<b>4001759 - 4001771 - 4002519</b>	2
<b>1000492*</b>	47	20	<b>4001772 - 4002518</b>	3

\*AISI 304 stainless steel blind baskets.

AISI 316 stainless steel wire dividers for autoclave baskets.

Part No.

**1001217** 4 divider compartments for basket part no. 1000495.

**1001218** 6 divider compartments for basket part no. 1000496.

**1001222** 7 divider compartments for basket part no. 1000780.

**Sterilization Drums.**

Código	Ø	Height cm	Suitable for autoclaves Part No.	Capacity Drums
<b>1002418</b>	25	18	<b>4001758 - 4001770 - 4002516</b>	3
<b>1002419</b>	35	27	<b>4001759 - 4001771 - 4002519</b>	2
<b>1002420</b>	28	18	<b>Not suitable for our autoclaves.</b>	

**Water distiller, specially for autoclaves “DEST-4”.**

Distillation capacity: 1.2 litres/hour. Internal capacity 4 litres.

Case made from thermo resistant plastic.

Internal stainless steel.

Cooling by forced air through a condenser. Without glass heater elements.

Easy to use. Safety cut out system if over temperature due to insufficient water.

Comes complete with a 4 litre three way dispenser plastic bottle with filter.

Note: for daily use it is recommended to change the filter every 3 months.



Part No.	Capacity Distillation litre/Hour	Conductivity µs/cm	Ø	Height cm	Power W	Weight Kg
<b>4001729</b>	1.2	5	29	39	750	3.5

**Spare:** Dispenser spout with filter. Part No. **4001730**

**Detergent “BIO-SEL”**

We recommend the use of BIO-Sel detergent for cleaning, it is also biodegradable and free from Phosphates.

Part No.

**1000025** Capacity 1 litre.

**1000026** Capacity 4 litres.

**1000852** Box of 12 x 1 litre flasks.

**1000853** Box of 4 x 4 litre flasks.