




without heating, analogical
[L] STIRRING



Mini magnetic stirrer 680

Code LBG001

The new mini magnetic stirrer Nahita-Blue is a basic laboratory instrument that presents a very dynamic design and an economic price. It is characterized for being compact and lightweight, which makes this stirrer an essential instrument for laboratories with limited space. It also comes with a vehicle power adapter, allowing it to be ideal for field work, where its use is required but without having a nearby power outlet.



Features

- It is a robust and durable product, as it is made of ABS plastic, allowing resisting minor splashing of almost any chemical.
- Simple to use, it has an analog control that allows agitation speed from 100 to 2000rpm.
- It is supplied with a support rod to use with probes for temperature, pH, ... and a stirring magnet, creating a fully functional and complete product for almost any need that the user may have.



Technical specifications

Code	LBG001
Model	680
Speed range	100-2000 rpm
Max. Volume	1 L
Plate dimensions	110 mm
Power	AC adapter 220V/12V
Dimensions (LxAxH)	120x120x50 mm
Weight	0,4 Kg

WITHOUT HEATING, ANALOGICAL
[L] STIRRING

Note: the manufacturer reserves the right to modify the characteristics and photographs of the products