



Cork supports

Code **BJQ001** - **BJQ002** - **BJQ003** - **BJQ004** - **BJQ005**

Available in several sizes to adapt to different flask volumes. Very stable supports used as stand for spherical flasks of different sizes.



Technical specifications

Code	BJQ001	BJQ002	BJQ003	BJQ004	BJQ005
Diameter	80 mm	100 mm	120 mm	140 mm	160 mm
d1	50 mm	65 mm	90 mm	110 mm	125 mm
d2	30 mm	45 mm	70 mm	90 mm	100 mm
H	30 mm	30 mm	30 mm	30 mm	30 mm
For flasks	50-100 mL	250 mL	500 mL	1000 mL	2000 mL

Rubber supports

Code **BJQ006**



With a maximum diameter of 90 mm it can hold flasks up to 500 mL.

Code	BJQ006
Diameter	90 mm
Height	50 mm

Polypropylene support

Code **BJQ007**



With a maximum diameter of 160 mm it can hold flasks up to 10 L

Code	BJQ007
Diameter	160 mm
Height	50 mm