



Stands




BOSSHEAD "P" TYPE

Nickel alloy bosshead.
For fixing perpendicular structures.
For 12 mm Ø bars.
Part No. **1000011**

ACCESSORY

Allen hex key. Suitable for "P" type bosshead.
Part No. **1000111**



J.P. Selecta

ORIGINAL BOSSHEAD "TREVINGE"

Made of chrome plated "Dural" alloy.
Due to its strength, suitable for multilevel and three dimension constructions.
Part No. **7008133** For Ø of 9 to 13 mm.
Part No. **7008253** For Ø of 19 to 25 mm.



RODS AND TUBES


12 mm Ø bars and rods. To be used with the "TREVINGE" and "P" type bossheads.

Dur-aluminium rods.

Part No.	Length
6000026	60 mm.
6000220	200 mm.
6000235	350 mm.
6000270	700 mm.
6002100	1000 mm.
6002150	1500 mm.
6002200	2000 mm.
6002300	3000 mm.

AISI 304 stainless steel tube.

Part No.	Length
6000016	60 mm.
6000120	200 mm.
6000135	350 mm.
6000170	700 mm.
6001100	1000 mm.
6001150	1500 mm.
6001200	2000 mm.
6001300 (rod)	3000 mm.



UNIVERSAL MOUNTING BASE

Bench top multiple rod support.
Suitable for mounting and supporting vertical and horizontal bars.
Part No. **1000050**



MOUNTING FOOT SUPPORT FOR RODS

For rod Ø of 12 mm.
Part No. **1000101** Rod position horizontal.
Part No. **1000049** Rod position vertical.



SIMULTANEOUS FOOT SUPPORT


For rod Ø of 12 mm.
Simultaneous positions, horizontal and vertical.
Part No. **1001102**



BENCH CLAMP "PASS"

Original Phywe. Clamps to the bench while holding from 4 to 12 mm diameter rod. Bench top width up to 40mm.
Part No. **1394855**


SUPPORT STAND KITS




700 high x 700 mm long.
 Made up of:
 8 rods 700 mm long.
 2 rods 350 mm long.
 16 "P" type bossheads.
 4 foot supports.
 1 Allen key.
 Part No. **6000090** Dur-aluminium
 Part No. **6000190** AISI 304 stainless steel



1000 high x 1000 mm long.
 Made up of:
 10 rods of 1000mm long.
 2 rods 700 mm long.
 3 rods 350 mm long.
 34 "P" type bossheads.
 6 foot supports.
 1 Allen key.
 Part No. **6000902** Dur-aluminium
 Part No. **6001902** AISI 304 stainless steel



700 high x 1000 mm long.
 Made up of:
 4 rods 1000 mm long.
 5 rods 750 mm long.
 2 rods 350 mm long.
 20 "P" type bossheads.
 4 foot supports.
 1 Allen key.
 Part No. **6000901** Dur-aluminium
 Part No. **6001901** AISI 304 stainless steel



1000 high x 2000 mm long.
 Made up of:
 5 rods 2000 mm long.
 7 rods 1000 mm long.
 8 rods 200 mm long.
 2 rods 350 mm long.
 45 "P" type bossheads.
 8 foot supports.
 1 Allen key.
 Part No. **6000903** Dur-aluminium
 Part No. **6001903** AISI 304 stainless steel




SUPPORT "STAND"
 Solid base, three vertical rods of 800 mm in height, two horizontal rods of 500 mm in length, six "P" type bossheads. Made of Dur-aluminium.
 Weight: 6.5 Kg.
 Part No. **6000482**


CENTRAL ROD RETORT STAND
 Vertical rod of Dur-aluminium of 10 Ø x 600mm high.
 Weight: 3 Kg.
 Base dimensions 315 x 200 mm.
 Part No. **6000046**



TRIPOD SUPPORT "PASS"
 Original Phywe.
 To support square and round rods from 6 to 14 mm dia.
 Weight: 2 Kg.
 Base with levelling screws.
 Part No. **6200255**



DOUBLE BASE SUPPORT T "MINI"
 Small foot support with screw clamp for removable rods from 8 to 12 mm Ø.
 Weight: 1.2 Kg.
 Part No. **6001103**



STANDARD RODS
 Made from AISI 340 stainless steel. Suitable for double "T MINI" and Tripod "PASS"
 Part No. **6001160** 8 Ø x 500 mm length
 Part No. **6000534** 10 Ø x 600 mm length
 Part No. **6000535** 12 Ø x 700 mm length



U SHAPED STAND

AISI 304 stainless steel of
20 Ø x 800 mm high.
Very stable.
Weight: 7 Kg.
Part No. **6000411**



DOUBLE BASE SUPPORT T "MAXI"

Rod of AISI 304 stainless steel of 20 Ø x 800 mm high.
Designed to support weights up to 20 Kg at 300 mm
height from the base. Useful for stirrers.
Weight: 7 Kg.
Part No. **6000410**



DOUBLE BASE SUPPORT T TELESCOPE

Stainless steel base AISI 304 of 32 mm Ø.
Telescopic elevation from 725 to 1025 mm.
Maximum support weight 20 Kg. Actual weight: 7.7 Kg
Compatible with Part No. 5057022 (See page 165).
Part No. **5057011**



TRIPOD SUPPORT STAND

AISI 304 stainless steel rod.

Part No.	Ø x Height	Weight Kg
6002278	8 x 500	1.4
6022710	10 x 600	1.8
6022712	12 x 700	2.5


Rods made of Dur-aluminium

Part No.	Ø x Height	Weight Kg
6000798	8 x 500	1.2
6000799	10 x 600	1.5
6007912	12 x 700	2



UNIVERSAL ARM SUPPORT

With holder for up to 4 sensors or elec-
trodes.
Very stable heavy base.
Ideal for working at varying heights.
Part No. **1001552**



RETORT STAND

Rods made of AISI 304 stainless steel

Part No.	Rod Ø x height	Plate size length x width	Weight Kg
6022614	10 x 600 mm.	215 x 130 mm.	3
6022612	12 x 700 mm.	315 x 200 mm.	3.6
6022613	12 x 800 mm.	315 x 200 mm.	3.8

Rods of Dur-aluminium


Part No.	Rod Ø x height	Plate size length x width	Weight Kg
6007714	10 x 600 mm.	215 x 130 mm.	2.7
6007712	12 x 700 mm.	315 x 200 mm.	3
6077812	12 x 800 mm.	315 x 200 mm.	3.4



TRIANGLE SUPPORTS


Made of 4 mm Ø AISI 304 stainless steel with electro-bright finish.

Part No.	Side	Part No.	Side
1002874	40 mm.	1002877	70 mm.
1002875	50 mm.	1002878	80 mm.
1002876	60 mm.	1028710	100 mm.



REFRACTORY TUBE AND WIRE TRIANGLE SUPPORTS

Part No.	Side
1000825	50 mm.
1000826	60 mm.
1000827	75 mm.



TRIPOD STAND "TRISOP"

Made of galvanised steel.

Part No.	Side	Total height
1000524	130 mm.	200 mm.



TRIPOD SUPPORTS WITH ROUND TOP

AISI 304 stainless steel

Part No.	Ø ext. ring	Total height
6002868	80 mm.	210 mm.
6028610	100 mm.	210 mm.
6028612	120 mm.	210 mm.
6028614	140 mm.	210 mm.

Nickel plated steel

Part No.	Ø ext. ring	Total height
1000818	80 mm.	210 mm.
1008110	100 mm.	210 mm.
1008112	120 mm.	210 mm.
1008114	140 mm.	210 mm.
1008116	160 mm.	210 mm.



TRIPOD STAND WITH RING TOP.

Galvanised steel.


Part No.	Ø ext. ring	Total height
7000808	80 mm.	210 mm.
7008010	100 mm.	210 mm.
7008012	120 mm.	210 mm.



TRIPOD STAND WITH RING TOP, ADJUSTABLE HEIGHT

Galvanised steel.

Part No.	Ø ext. ring	Adjustable height
7022810	100 mm.	170 to 250 mm.
7022812	140 mm.	170 to 250 mm.
7022814	180 mm.	170 to 250 mm.



1. DECANTING FLASK SUPPORT

Made of polypropylene. With support flasks up to 12 cm Ø.

Part No. **1001346**

2. DECANTING FLASK SUPPORT WITH BOSSHEAD

Made of polypropylene. With support flasks from 30 to 150 mm Ø. Complete with reduction ring.

Part No. **1001347**