

▶ Beaker PP low form

- 1 | Made of polypropylene.
- 2 | With spout to prevent spillage.
- 3 | With permanent embossed graduation.



Code	Capacity	Graduation	Resolution	Height	Diameter
CBH001	50 ml	10 ml	5 ml	60 mm	42 mm
CBH002	100 ml	20 ml	5 ml	70 mm	54 mm
CBH003	250 ml	50 ml	10 ml	100 mm	71 mm
CBH004	500 ml	100 ml	10 ml	120 mm	92 mm
CBH005	1000 ml	100 ml	20 ml	145 mm	110 mm