




► Groove spoon spatulas



DIMENSIONS - groove spoon spatulas




Code	L (mm)	AxB (mm)	CxD (mm)
FKB015	120	15x13	15x10
FKB016	140	20x15	17x11
FKB017	160	25x18	17x12
FKB018	180	30x20	20x13

► Spoon spatulas



DIMENSIONS - spoon spatulas



Code	L (mm)	AxB (mm)
FKP003	125	20x13
FKP004	160	25x15
FKP005	180	25x18
FKP006	220	30x18

► Filled spatula



DIMENSIONS - filled spatula




Code	AxB (mm)	CxD (mm)	L (mm)
FKP001	25x11	50x9	215

► Double spoon spatulas



DIMENSIONS - double spoon spatulas



Code	L (mm)	AxB (mm)	CxD (mm)
FKR001	160	25x18	15x11
FKR002	200	30x20	25x15



► Spoon-flat spatulas



DIMENSIONS - spoon-flat spatulas




Code	L (mm)	AxB (mm)	CxD (mm)
FKB005	150	15x20	25x12
FKB006	200	30x20	32x15
FKB007	230	30x20	30x15

► Spoon-flat spatulas, reinforced



DIMENSIONS - spoon-flat spatulas, reinforced




Code	L (mm)	AxB (mm)	CxD (mm)
FKB008	100 mm	20x15 mm	22x13 mm
FKB009	125 mm	20x14 mm	18x14 mm
FKB010	150 mm	25x19 mm	25x19 mm
FKB011	210 mm	30x25 mm	30x25 mm
FKB012	250 mm	28x21 mm	35x20 mm

► Flat-curved spoon spatula



DIMENSIONS - flat-curved spoon spatula



Code	L (mm)	AxB (mm)	CxD (mm)
FKB020	195	20x13	30x20

► Analysis spatula with spoon



DIMENSIONS -analysis spatula with spoon




Code	L (mm)	AxB (mm)	CxD (mm)
FKB021	224	40x19	48x10



► Trough spatula, groove-spoon



DIMENSIONS - trough spatula, groove-spoon



Code	L (mm)	AxB (mm)	CxD (mm)
FKB022	200	30x24	90x17

► Special spoon spatula



Code	With (mm)	Lenght (mm)
FKB042	7	60

► Microspatula spoon-flat



DIMENSIONS - microspatula spoon-flat




Code	L (mm)	AxB (mm)	CxD (mm)
FKD001	190	35x9	50x9
FKD002	240	35x16	52x9
FKD003	150	12x8	54x7

► Spatula flat-flat reinforced



DIMENSIONS - spatula flat-flat reinforced




Ref. / Code	AxB (mm)	L (mm)
FKL001	35x18	150
FKL002	45x18	200



► Chattaway spatulas



DIMENSIONS - chattaway spatulas



Code	L (mm)	AxB (mm)
FKN001	100	25x8
FKN002	120	25x9
FKN003	150	25x9
FKN004	170	30x9
FKN006	140	35x7

► Microspatula flat-groove




Code	Spatula AxB (mm)	L (mm)
FKN007	30x5	150
FKN008	30x5	200

► Spatulas with wooden handle



DIMENSIONS - spatulas with wooden handle



Code	AxB (mm)	L (mm)
FKG006	100x20	215
FKG007	150x25	270
FKG009	200x25	315