

Papel para lentes | **Lens tissues**

Marca|**Trademark**: Ura Technic



Referencia / **Code** : 89000000

Las lentes y otras superficies ópticas de vidrio, plástico o cuarzo pueden ser fácilmente rayadas si no se limpian adecuadamente. El papel para lentes Zuzi, ligero y sin pelusas, proporciona la solución ideal para la eliminación de polvo, humedad o grasa presente en las superficies ópticas de microscopios, lupas, cámaras fotográficas y otros equipos delicados de laboratorio. Se recomienda que siempre se aplique el papel para lentes humedecido con una solución limpiadora para evitar dañar la superficie óptica.

Cada pliego de papel para lentes Zuzi contiene alrededor de 50 unidades de 110x85 mm de tamaño.

Lenses and other optical surfaces made of glass, plastic or quartz can be easily scratched and damaged if they are not cleaned properly. Zuzi lens tissues, light and lint free, are the ideal solution to wipe away dust, moisture and grease from optical surfaces of microscopes, magnifying glasses, photographic cameras and other delicate instruments of laboratory. It is recommended to apply the lens tissue previously moistened with a cleaning solution to avoid damages on the optical surface.

Each wallet of Zuzi lens tissues contains approximately 50 sheets of 110x85 mm.

