



FUNNEL WITH HANDLE

Made of AISI304 stainless steel.

Part No. 1000906	100 mm Ø x 100 mm high x 10 mm exit.
Part No. 1000907	120 mm Ø x 112 mm high x 10 mm exit.
Part No. 1000908	150 mm Ø x 159 mm high x 20 mm exit.
Part No. 1000909	200 mm Ø x 183 mm high x 20 mm exit.



SPOON WITH HANDLE

Made of AISI304 stainless steel.

Part No. 1001594	Capacity 100 ml.	Total length 240 mm. Ø 83 mm
Part No. 1001595	Capacity 200 ml.	Total length 280 mm. Ø 92 mm
Part No. 1001596	Capacity 300 ml.	Total length 360 mm. Ø 102 mm



FLAT BOTTOM BASINS

Stainless steel AISI 304	Nickel Ni	Ø	Height	Volume
Part No. 1009160	Part No. 1019160	60 mm.	25 mm.	40 ml.
1009170	1019170	70 mm.	25 mm.	60 ml.
1009180	1019180	80 mm.	30 mm.	90 ml.
1091100	1191100	100 mm.	35 mm.	180 ml.



EVAPORATION BASINS, ROUND BOTTOM

Stainless steel AISI 304	Nickel Ni	Ø	Height	Volume
Part No. 1009260	Part No. 1019260	60 mm.	25 mm.	35 ml.
1009270	1019270	70 mm.	25 mm.	50 ml.
1009280	-	80 mm.	30 mm.	75 ml.
1092100	-	100 mm.	35 mm.	125 ml.



PETRI DISH WITH LID

Stainless steel AISI 304	Nickel Ni	Ø	Height	Volume
Part No. 1002126	Part No. 1012126	60 mm.	20 mm.	50 ml.
1002128	1012128	80 mm.	20 mm.	90 ml.
1021210	-	100 mm.	25 mm.	150 ml.
1021212	-	120 mm.	25 mm.	260 ml.
1021215	-	150 mm.	30 mm.	550 ml.



BUCKETS WITH HANDLES

Without graduation. Made of AISI304 stainless steel.

Part No.	Height	Ø	Capacity
1001599	270 mm	270 mm	10 litres
1001600	320 mm	350 mm	19 litres

BUCKET LIDS

Part No. 1001601	Suitable for buckets 10 litres Part No. 1001599
Part No. 1001602	Suitable for buckets 19 Litres Part No. 1001600