

Pipettors Compatibility Testing

Pipettors Compatibility Testing		10 µl			20 µl		100 µl		200 µl		1.000 µl										
		Non filter tips			Filter tips	Filter tips		Filter tips	Non filter tips		Filter tips	Non filter tips			Filter tips						
		MJB032	MJB034	MJB035	MJB033	MJB036	MJD010	MJD011	MJD012	MJC005	MJC006	MJJ040	MJJ042	MJJ043	MJJ041	MJJ044	MJN038	MJN040	MJN041	MJN039	MJN042
BioHit																					
	BM-3 Biohit mLine manual pipettor, 1 Ch. 0.1-3ul	●	●	●																	
	BM-10 Biohit mLine manual pipettor, 1 Ch. 1-10ul	●	●	●																	
	BM-20 Biohit mLine manual pipettor, 1 Ch. 2-20ul						●	●	●			●	●	●							
	BM-100 Biohit mLine manual pipettor, 1 Ch. 10-100ul									●	●										
	BM-200 Biohit mLine manual pipettor, 1 Ch. 20-200ul											●	●	●							
	BM-1000 Biohit mLine manual pipettor, 1 Ch. 100-1000ul																●	●	●		●
	BU-20 Proline Plus 1 channel pipettor 2-20ul						●	●	●			●	●	●							
	BU-100 Proline Plus 1 channel pipettor 10-100ul									●	●										
	BU-200 Proline Plus 1 channel pipettor 20-200ul											●	●	●							
	BU-1000 Proline Plus 1 channel pipettor 100-1000ul																●	●	●		●
	BH-20 Biohit Proline Manual Pipettor, 1 ch., 2-20ul						●	●	●			●	●	●							
	BH-50 Biohit Proline Manual Pipettor, 1 ch., 5-50ul											●	●	●							
	BH-100 Biohit Proline Manual Pipettor, 1 ch., 10-100ul									●	●	●	●	●							
	BH-200 Biohit Proline Manual Pipettor, 1 ch., 20-200ul											●	●	●		●	●				
	BH-1000 Biohit Proline Manual Pipettor, 1 ch., 100-1000ul																●	●	●		●
	BH8-50 Biohit Proline Manual Pipettor. 8 ch., 5-50ul in 0.5ul											●	●	●				●	●	●	●
	BL-120A 1 ch. 5-120ul w/Adapter											●	●	●							
	BP-100 Single channel, 5-100ul in 1ul increments									●	●	●	●	●							
	BP-500 Single channel, 10-500ul in 5ul increments											●	●	●							
	BP-1000 Single channel, 50-1000ul in 10ul increments																●	●	●	●	●
Capp	BP12-100 Twelve channel, 5-100ul in 1ul increments											●	●	●		●	●				
	CU-2 CappAero 1 ch. pipette, 0.2-2 ul Standard Gray	●	●	●																	
	CU-10 CappAero 1 ch. pipette, 1-10 ul Standard Gray	●	●	●																	
	CU-20 CappAero 1 ch. pipette, 2-20 ul Standard Gray						●	●	●			●	●	●							
	CU-50 CappAero 1 ch. pipette, 5-50 ul, Ecoline, Red											●	●	●							
	CU-100 CappAero 1 ch. pipette, 10-100 ul, Ecoline, Orange									●	●										
	CU-200 CappAero 1 ch. pipette, 20-200 ul, Ecoline, Yellow											●	●	●		●	●				
Eppendorf	CU-1000 CappAero 1 ch. pipette, 100-1000 ul, Ecoline, Blue																●	●	●		●
	E2000-10 1 ch. from 1-10ul Air-Displacement Micropipette																				
	E2000-20 1 ch. from 2-20ul Air-Displacement Micropipette																				
	E2000-100 1 ch. from 10-100ul Air-Displacement Micropipette								●	●		●	●	●							
	E2000-200 1 ch. from 50-200ul Air-Displacement Micropipette											●	●	●		●	●				
	E2000-1000 1 ch. from 100-1000ul Air-Displacement Micropipette																●	●	●		●
	E2100-2.5 1 ch. from 0.1-2.5ul Air-Displacement Micropipette																				
	E2100-10 1 ch. from 1-10ul Air-Displacement Micropipette																				
	E2100-20 1 ch. from 2-20ul Air-Displacement Micropipette						●	●	●			●	●	●							
	E2100-100 1 ch. from 10-100ul Air-Displacement Micropipette									●	●										
	E2100-200 1 ch. from 20-200ul Air-Displacement Micropipette											●	●	●		●	●				
	E2100-1000 1 ch. from 100-1000ul Air-Displacement Micropipette																●	●	●		●
	ES-1000 Research Plus																●	●	●		●
	ES-100 Research Plus 10-100ul									●	●	●	●	●							
	ES-200 Research Plus 20-200ul																				
	EE-100 1 ch. from 10-100ul Electronic Multichannel Pipette									●	●	●	●	●							
	EE-1000 1 ch. from 100-1000ul Electronic Pipette																●	●	●		●
	EE-10																				
Fisher Scientific	EE12-10																				
	FFI-50 FisherBrand Finnpipette I 5-50 ul											●	●	●							
	FFI-200 FisherBrand Finnpipette I 20-200 ul											●	●	●		●	●				
Gilson	FFI-1000 FisherBrand Finnpipette I 100-1000 ul																●	●	●		●
	P-10 Pipetman 1 ch. Air-Displacement Micropipette , 1-10 ul	●	●	●	●	●															
	P-20 Pipetman 1 ch. Air-Displacement Micropipette , 2-20 ul						●	●	●			●	●	●							
	P-20 Pipetman multi ch. Air-Displacement Micropipette																				
	P-100 Pipetman 1 ch. Air-Displacement Micropipette , 10-100 ul									●	●	●	●	●							
	P-200 Pipetman 1 ch. Air-Displacement Micropipette , 50-200 ul											●	●	●		●	●				
	P-200 Pipetman Multichannel Air-Displacement Micropipette																				
Labnet/HTL	P-1000 Pipetman 1 ch. Air-Displacement Micropipette , 100-1000 ul																				
	LN-2 Labnet Labpette Discovery 0.2 to 2ul autoclavable	●	●	●	●	●															
	LN-20 Labnet Labpette Discovery 2 to 20ul autoclavable						●	●	●			●	●	●							
	LN-100 Labnet Labpette Discovery 10 to 100ul autoclavable									●	●										
Nichiryo	LN-200 Labnet Labpette Discovery 20 to 200ul autoclavable											●	●	●		●	●				
	LN-1000 Labnet Labpette Discovery 100 to 1000ul autoclavable																●	●	●		●
	NPX-2 Nichipet EX Autoclavable/UV stable digital 0.2-2ul	●	●	●	●	●															
	NPX-10 Nichipet EX Autoclavable/UV stable digital 1-10ul	●	●	●	●	●															
	NPX-20 Nichipet EX Autoclavable/UV stable digital 2-20ul						●	●	●			●	●	●							
	NPX-100 Nichipet EX Autoclavable/UV stable digital 10-100ul									●	●	●	●	●							
	NPX-200 Nichipet EX Autoclavable/UV stable digital 20-200ul											●	●	●							
Rainin	NPX-1000 Nichipet EX Autoclavable/UV stable digital 100-1000ul																				
	ED1-200R Rainin EDP 1 ch.Electronic Pipette , 20 - 200 ul											●	●	●		●	●				
	ED2-25R Rainin EDP 2 1 ch. Electronic Pipette , 2.5 - 25 ul											●	●	●							
	ED2-100R Rainin EDP 2 1 ch. Electronic Pipette , 10 - 100 ul									●	●	●	●	●							
	ED2-1000R Rainin EDP 2 1 ch. Electronic Pipette , 100 - 1000 ul																●	●	●		●
	ED-100R Rainin EDP Plus ch. Electronic Pipette , 10 - 100 ul									●	●	●	●	●							
Socorex	ED-1000R Rainin EDP Plus ch. Electronic Pipette , 100 - 1000 ul																				
	S825-20 Acura® 825 Single-Ch. Variable Volume 2-20						●	●	●			●	●	●							
	S825-50 Acura® 825 Single-Ch. Variable Volume 5-50											●	●	●							
	S825-100 Acura® 825 Single-Ch.Variable Volume10-100									●	●	●	●	●							
	S825-200 Acura® 825 Single-Ch. Variable Volume 20-200											●	●	●		●	●				
	S825-1000 Acura® 825 Single-Ch. Variable Volume 100-1000																●	●	●		●