



Glass syringe Luer

- Made of neutral glass. Plunger tip marked with amber ring
- Interchangeable. Individually packaged .
- Steam sterilizable by autoclave. It is recommended to do it at 134°C for 10 minutes, separating the piston from the cylinder.
- Boquilla Luer de vidrio, Luer metal o luer
NOT FOR MEDICAL USE

Technical especificacions

Capacity	Luer Lock	Luer metal	Luer glass
1 mL	MMB002	MMB011	MMB043
2 mL	MMB003	MMB032	MMB044
3 mL	MMB004	MMB033	MMB045
5 mL	MMB005	MMB034	MMB046
10 mL	MMB006	MMB013	-
20 mL	MMB007	MMB012	-
30 mL	MMB008	MMB014	MMB049
50 mL	MMB009	MMB035	MMB050
100 mL	MMB010	-	-