

Médula espinal / Spinal cord

Aparecen representadas las sustancia gris (astas anterior y posterior, comisura gris y conducto epéndimo), sustancia blanca (cordones anterior, lateral y posterior, comisura anterior, surcos ventral y dorsal) y las raíces motora y sensitiva. Se suministra con base

Representing grey matter (anterior and posterior horns, grey commissure and central canal), white matter (anterior, lateral and posterior funiculus, anterior commissure, anterior median fissure and posterior median sulcus) as well as anterior and posterior roots. Supplied with a baseboard.




Referencia / Code	48000170
Partes / Parts	1
Aumentos / Magnification	2x

Cerebro / Brain

Dividido por la mitad, el modelo 48000190 se puede desmontar en lo lóbulos frontal y parietal, lóbulo temporal con occipital, mitad del tallo cerebral y mitad del cerebelo. Se suministra con base

Divided in two parts, model 48000190 can also be disassembled into frontal and parietal lobes, temporal and occipital lobes, half the brainstem and half the cerebellum. Supplied with a baseboard.



Referencia / Code	48000190	48000200
Aumentos / Magnification	1x	1x
Partes / Parts	9	2

Nota: el fabricante se reserva el derecho de modificar las características y fotografías de los productos.
Note: Manufacturer reserves the right to modify features and photographs of the goods.

