



### ▶ Narrow bottle, PE



- 1 | Round bottle with narrow neck without cap
- 2 | Made in polyethylene (PE).
- 3 | Natural colour

| Code   | GL | Cap.   | Diameter | Ø neck | Height |
|--------|----|--------|----------|--------|--------|
| CGJ047 | 18 | 50mL   | 40 mm    | 13 mm  | 84 mm  |
| CGJ048 | 18 | 100mL  | 45 mm    | 13 mm  | 102 mm |
| CGJ056 | 25 | 250mL  | 60 mm    | 18 mm  | 138 mm |
| CGJ049 | 25 | 500mL  | 75 mm    | 18 mm  | 169 mm |
| CGJ050 | 28 | 1000mL | 95 mm    | 21 mm  | 212 mm |
| CGJ051 | 28 | 2000mL | 120 mm   | 21 mm  | 260 mm |

### ▶ Screw caps, PE



- 1 | Provided with holes for security seals. Made of PP.

| Code   | GL |
|--------|----|
| CGN011 | 18 |
| CGN012 | 25 |
| CGN013 | 28 |

### ▶ Wash bottle screw caps



- 1 | For standardized narrow neck bottles.

| Code   | DIN |
|--------|-----|
| CKL004 | 18  |
| CKL005 | 25  |
| CKL006 | 28  |

### ▶ Short spout dropping



- 1 | For standardized narrow neck bottles.

| Code   | Spout type | GL |
|--------|------------|----|
| CKL007 | Short      | 18 |
| CKL008 | Short      | 25 |
| CKL009 | Long       | 18 |
| CKL010 | Long       | 25 |
| CKL011 | Long       | 28 |