


► Frascos PP boca ancha

- 1| Cuello ancho, PP.
- 2| Autoclavables a 121°C.
- 3| Con tapa.
- 4| Color natural.



Referencia	Capacidad	Diámetro	Diámetro boca	Altura	GL
CGD071	50 ml	40 mm	24 mm	77 mm	32
CGD072	100 ml	48 mm	24 mm	92 mm	32
CGD073	250 ml	64 mm	31 mm	123 mm	40
CGD074	500 ml	75 mm	38.5 mm	153 mm	50
CGD075	1000 ml	94 mm	50 mm	202 mm	65
CGD076	2000 ml	120 mm	50 mm	243 mm	65