



Tube roller, RSLAB-10



Features

It allows shaking several types of tubes with different liquid, viscous samples or solids suspension that require minimum airing.

- 1 | 6 rollers.
- 2 Rocking and rolling action for complete mixing.

- 3 | Analogical regulation speed.
- 4 Continuous operation mode.
- 5 | Protection class according to DIN EN 60529: IP 21

Technical features

Code	LJD002
Model	RSLAB-10
Shaking movement	Rotary and rolling
Speed range	0-70 rpm
Max. shaking weight	5 Kg
Run type	Time/continuous
Dimensions (LxWxH)	530x260x120 mm
Weight	5.1 Kg

MULTIDIMENSIONAL SHAKERS [L] STIRRING

Note: manufacturer reserves the right to modify features and photographies of the goods.