




Spatulas, scoops and vessels.

STAINLESS STEEL AND NICKEL EVAPORATING BASINS.
STAINLESS STEEL, NICKEL AND ZIRCONIUM CRUCIBLES.



MICRO SPATULA WITH FLAT ENDS

Made of AISI 304 stainless steel.


Part No.	Length	Width
1027680	80 mm.	4 mm.
1276120	120 mm.	4 mm.
1276150	150 mm.	4 mm.
1276210	210 mm.	4 mm.



MICRO SPATULA WITH FLAT AND SPOON ENDS

Made of AISI 304 stainless steel.

Part No.	Length	Width
1277100	100 mm.	4 mm.
1277150	150 mm.	4 mm.
1277210	210 mm.	4 mm.



MICRO-SPATULA CURVED ENDS

Made of stainless steel AISI 304.

Part No. **1001581** Length 150 mm.



SPATULA FLAT AND GROOVED ENDS

Made of Stainless steel AISI 304.

Part No. **1000855** Total length 180 mm. Flat width 9 mm. Grooved width 8 mm.



SPATULA WITH CURVED ENDS

Made of stainless steel AISI 304.

Part No.	Length	Width
1010313	130 mm.	11 mm.
1001518	150 mm.	11 mm.
1010318	180 mm.	11 mm.
1010321	210 mm.	11 mm.



SPATULA WITH SPOON AND FLAT ENDS

Part No. Length Width


1096120	120 mm.	18 mm.	stainless steel AISI 304.
1096150	150 mm.	20 mm.	stainless steel AISI 304.
1196120	120 mm.	18 mm.	Pure nickel.
1196150	150 mm.	20 mm.	Pure nickel.



CHISEL AND GROOVED SPOON END SPATULA

Part No. Length Width


6021312	120 mm.	18 mm.	stainless steel AISI 304.
6021315	150 mm.	20 mm.	stainless steel AISI 304.
6121312	120 mm.	18 mm.	Pure nickel.
6121315	150 mm.	20 mm.	Pure nickel.



FLAT SPATULA WITH SPOON

Made of Stainless steel AISI 304.

Part No. **1001583** Length 180 mm.



SPATULA DOUBLE SPOON ENDS

Made of Stainless steel AISI 304.

Part No. **1001582** Length 150 mm.




DOUBLE SPATULA SEMI-MICRO

Made of Stainless steel AISI 304. Total length 140 mm

Rectangular end: 35 long x 7 mm wide.

Oval end: 35 long x 7 mm wide.

Part No. **1278140**



SPATULA WITH POCKET CLIP

Part No.

1. **1001585** Micro-spatula flat spoon end. Length 150 mm.
2. **1001586** Micro-spatula curved end. Length 150 mm.
3. **1001587** Spatula flat with grooved end. Length 150 mm.
4. **1001588** Spatula double end flat. Length 150 mm.
5. **1001589** Spatula flat and spoon ends. Length 150 mm.



Specially designed for lab coats with upper pocket.



FLAT AND CURVED END SPATULA

Made of AISI 304, stainless steel.


Part No.	Length	Width
6099120	120 mm.	20 mm.
6099150	150 mm.	24 mm.
6099180	180 mm.	24 mm.
6099210	210 mm.	24 mm.



SPATULA DOUBLE CURVED ENDS

Made of AISI 304, stainless steel.

Part No.	Length	Width
6100120	120 mm.	15 mm.
6100150	150 mm.	25 mm.
6100180	180 mm.	25 mm.
6100210	210 mm.	25 mm.



GROOVED SPATULA

Grooved spoon, made of AISI 304 stainless steel.

Total length: 190 mm.

Designed to remove residues from flasks and bottles due to the flat sharp end and curved spoon.

Part No. **1096190**



SPATULA SPOON FLAT-CONCAVE

Made Stainless steel AISI 304.

Part No. **1001584** Length 250 mm.



FLAT SPATULA WITH SPOON

Analytical type, made of AISI 304 Stainless steel. Total length: 210 mm.

Rectangular end: 40 long x 9 mm wide.

Spoon end: 30 long x 18 mm wide.


Part No. **1095210**



PASTE KNIFE WITH WOODEN HANDLE

Blade made from tempered Stainless steel.

Part No.	Blade length	Blade width	Part No.	Blade length	Blade width
11010280	80 mm.	14 mm.	1102160	160 mm.	24 mm.
1102100	100 mm.	16 mm.	1102200	200 mm.	29 mm.
1102120	120 mm.	18 mm.	1102240	240 mm.	34 mm.
1102140	140 mm.	22 mm.	1102280	280 mm.	34 mm.



FLEXIBLE SPATULA

With stainless steel AISI 304 blade and handle.

Part No.	Total length	Blade length	Blade width
1000856	160 mm.	52 mm.	11 mm.
1000857	160 mm.	52 mm.	15 mm.
1000858	170 mm.	60 mm.	20 mm.



MULTI PURPOSE KNIFE

With a black bakelite® handle. Designed for comfort and safety
Made of tempered stainless steel.
Total length: 230 mm.
Blade length: 95 mm.
Part No. **1000862**



SPATULA POCKET POUCH

Part No. **1001530**

Contents:

1 Flexible spatula, Part No. 1000856 . 1 Semi micro spatula, Part No. 1278140
1 Micro spatula flat spoon, Part No. 1277150. 1 Dissection scissors, Part No. 1000807
1 Forceps blunt end, Part No. 1000580



SODIUM SPOON

With brass basket. Ø 25 mm. Total length: 190 mm.
Part No. **1000805**




POWDER SCOOPS

Length: 100 mm. Width: 26 mm.

Part No. **6000278** Stainless steel AISI 304.

Part No. **6001278** Pure nickel.



SCOOPS WITH HANDLE

Made of stainless steel AISI 304.

Part No. **1001592** Medium Total length 320 mm.

Part No. **1001593** Big Total length 360 mm.



VIBRATING SPATULA AND BALANCE WEIGHT FORCEPS

Vibrating spatula.

With manual vibrating device. Made of AISI 304 stainless steel.

Part No. **1000413** Curve end. Total length 176 mm. Width 12 mm.

Part No. **1000412** Grooved end. Total length 176 mm. Width 12 mm.

Balance weight forceps.

Made of AISI 304 stainless, electrobright finish. Total length 93 mm.

Part No. **1000217**

Forceps for balance weights.

Made of AISI 304 stainless, electrobright finish. Total length 100 mm.

Part No. **1000216**



FUSION SPOONS

Made of AISI 304 stainless, wooden handle. Total length: 350 mm.

Flat base. Volume 90 ml.

Part No. **6000093**

Round base. Volume 75 ml.

Part No. **6000094**



SCOOPS WITH HANDLE

Made of stainless steel AISI 304.

Part No.	Total length	Usable length	Width	Capacity
1000859	180 mm.	100 mm.	50 mm.	120 ml.
1000860	260 mm.	150 mm.	75 mm.	450 ml.
1000861	310 mm.	200 mm.	100 mm.	1500 ml.
1001354	330 mm.	210 mm.	150 mm.	2000 ml.